


Scientific CMOS

Zyla and Neo sCMOS Cameras Widen Your Expectations



Scientific CMOS (sCMOS)

Scientific CMOS (sCMOS) is a breakthrough technology that offers an advanced set of performance features that render it ideal to high fidelity, quantitative scientific measurement. Scientific CMOS can be considered unique in its ability to simultaneously deliver on many key performance parameters, overcoming the 'mutual exclusivity' associated with current scientific imaging technology standards, and eradicating the performance drawbacks traditionally associated with CMOS imagers.

The multi-megapixel sensors offer a large field of view and high resolution, without compromising read noise, dynamic range or frame rate. Rolling and Global (Snapshot) shutter readout ensure maximum application flexibility.

Read noise is exceptional, even when compared to the highest performance 'slow-scan' CCDs. The fact that an sCMOS device can achieve 1 electron rms read noise while reading out up to 5.5 megapixels at 100 fps through the 10-tap cameralink interface, or 40 fps through the USB 3.0 interface renders it truly extraordinary in the market. By way of comparison, the lowest noise Interline CCD, reading out only 1.4 megapixels at ~16 fps would do so with ~10 electrons read noise.

The low noise readout is complemented by up to 33,000:1 dynamic range. Usually, for CCDs or EMCCDs to reach their highest dynamic range values, there needs to be a significant compromise in readout speed, yet sCMOS can achieve this value while delivering high frame rates. The unique dual amplifier architecture of sCMOS allows for high dynamic range by offering a large well depth, despite the relatively small 6.5 µm pixel size, alongside lowest noise. A 1.4 megapixel Interline CCD with similarly small pixels achieves only ~1,800:1 dynamic range at 16 fps.

	P	

Neo cameras will literally allow one to see cells in a new light with ultra sensitive imaging at speeds never achieved before - as we have seen in our tests of vesicle trafficking. These scientific CMOS cameras are not a small step, but a quantum leap, that will open up new possibilities of what can be studied in fast cellular processes, rapid screening, and super-resolution imaging.

Derek Toomre, PhD., Associate Professor, Department of Cell Biology, Yale University School of Medicine, USA

	sCMOS - Zyla 4.2 PLUS
	4.2 Megapixel
	6.5 µm
	0.9e- @ 30 fps / 1.1 e- @100 fps
(max)	Sustained: 100 fps full frame
ncy (max)	82%
	33,000:1 (@ 30 fps)
bise	None

Parameter

Sensor Forma

Full Frame Rate

Quantum Efficie

Dynamic Range

Multiplicative N

Pixel Size <u>Re</u>ad Noise



Key sCMOS Features

Extremely low noise

Rapid frame rates

Wide dynamic range

High resolution

Large field of view

High Quantum Efficiency (QE)

Rolling and Global (Snapshot) exposure modes



Neo 5.5 sCMOS

Driving optimal imaging performance



Andor's Neo sCMOS vacuum cooled camera platform has been engineered from the ground up, specifically to realize the absolute highest sensitivity from this exciting new sensor technology.

Neo 5.5 offers an exceptionally low dark current and read noise floor detection limit, maintained even under longer acquisition times, alongside a wide dynamic range of 30,000:1. Speeds of 30 fps (full frame) can be maintained over extended kinetic series acquisitions, with 100 fps achievable in burst exposure simultaneously. mode.

Neo 5.5 offers an advanced set of unique performance features and innovations, including deep vacuum TE cooling to -40°C, extensive 'on-head' FPGA data processing capability, a 4 GB memory buffer and a Data Flow Monitor. Andor's UltraVacTM vacuum process has been implemented to offer not only the necessary deep cooling capability, but also complete protection of the sensor.

been conceptualized to drive the best possible performance, image quality and longevity from sCMOS technology.

Neo 5.5 offers both Rolling and true Global (also known as 'Snapshot') shutter exposure mechanisms. Snapshot mode provides an exposure sequence that is analogous to that of an Interline CCD, whereby all pixels begin the exposure simultaneously and end the



andor.com/neo

Key Specifications

Active pixels	2560 x 2160
Pixel size (w x h; µm)	6.5 x 6.5
Image area (mm)	16.6 x 14
Active area pixel well dep	oth (e) 30,000
Max readout rate (MHz)	560 (280 x 2 outputs)
Frame rates (fps) S	Sustained: 30 (full frame) Burst: 100 (full frame)
Read noise (e ⁻)	1 @ 200 MHz 1.3 @ 560 MHz
Sensor cooling	-40°C
Correct Coolining	40.0

Features and Benefits

TE cooling to -40°C	Minimization of dark current of hot pixel blemishes mear
Rolling and true Global (Snapshot) shutter	Maximum exposure and rea capability.
1 e- read noise	Offers lower detection limit t
5.5 megapixel sensor format and 6.5 µm pixels	Delivers extremely sharp res pathology, high content scre
Dark Noise Suppression (DNS) technology	Extremely competitive low c advantage across range of
Rapid frame rates	>30 fps over extended kine
UltraVac™	Sustained vacuum integrity
Dual-Gain amplifiers	Maximum well depth and lo
GPU Express	Simplify and optimize data t accelerated GPU processin
4 GB on-head image buffer	Enables bursts of 100 fps @ speed, avoiding prohibitively
Sub-microsecond inter-frame gap	Global Shutter offers down
Extensive FPGA on-head data processing	Essential to ensure best ima
Hardware timestamp	FPGA generated timestamp
Dynamic baseline clamp	Essential to ensure quantita kinetic series.
Spurious noise filter	Realtime FPGA filter that ide
Data flow monitor	Innovatively manage acquis
Comprehensive trigger modes and I/O	Communication and synchr

Key Applications
Neuroimaging
Cell Motility
Veolcimetry
Super-resolution
TIRF Microscopy
FRET / FRAP
Light sheet microscopy
Vesicle Trafficking
Ion Signaling
Astronomy

Neo 5.5 QE curve

	70	
	70	
	60	
	50	
(%)	40	
Quantum efficiency (%)		
antum €	30	
ğ	20	
	10	
	0	00
	3	00



nt to maintain low noise advantage under all exposure conditions. Minimization aning more useful pixels. Fan-off mode for vibration sensitive set-ups.

eadout flexibility across all applications. Snapshot for 'Interline CCD' exposure

than any CCD.

esolution over a 22 mm diagonal field of view: ideal for cell microscopy, digital reening and astronomy.

dark current of 0.015 e⁻/pix/sec with fan cooling. Maintains low noise f exposure conditions.

etic series. Burst to memory at 100 fps full frame.

and unequalled cooling with 5 year warranty; complete sensor protection.

lowest noise simultaneously, affording extended dynamic range of 30,000:1.

transfers from camera to Graphical Processing Unit (GPU) card to facilitate ng as part of the acquisition pipeline.

@ full dynamic range. Capture extended kinetic series faster than PC write ly expensive PCs.

n to 100 ns inter-frame gap, ideal for PIV applications.

nage quality and quantitative fidelity from sCMOS technology.


np with 25 ns accuracy.

ative accuracy across the image area and between successive images of a

dentifies and compensates for spurious high noise pixels.

sition capture rates vs data bandwidth limitations.

nronization within intricate experimental set-ups.



Zyla 5.5 sCMOS

High resolution camera for research and OEM



Andor's Zyla 5.5 sCMOS camera offers high speed, high sensitivity imaging in a remarkably light and compact design. Both Rolling and true Global shutter modes offer extensive application flexibility. Global shutter is ideally suited to fast multi-dimensional microscopy, offering tight synchronization to 'moving' peripheral devices such as z-stage or light source.

Zyla is ideally suited to many cutting edge applications that push the boundaries of speed, offering sustained frame rate performance of up to 100 fps (faster with ROI).

A highly cost-effective USB 3.0 version is also available, offering an unparalleled 40 fps (full frame) and 0.9 e- rms read noise, representing an ideal low light 'workhorse' camera solution for both microscopy and physical science applications, in either research or OEM environments. Rolling and Global (Snapshot) shutter readout ensures maximum application flexibility. Global shutter in particular provides an important 'freeze frame' exposure mechanism that emulates that of an Interline CCD, overcoming the transient readout nature of Rolling shutter mode.

Zyla sCMOS for OEM

The light and compact form factor coupled with design and mounting adaptability, board level or private labelling options, and unparalleled engineering support, renders the Zyla highly suited to OEM integration.

Please call Andor to discuss how Zyla can be made to work for you.

Key Specifications

Active pixels	2560 x 2160
Pixel size (w x h; µm)	6.5 x 6.5
Image area (mm)	16.6 x 14
Active area pixel well dep	oth (e) 30,000
Max readout rate (MHz)	560 (280 x 2 outputs)
Frame rates (fps)	Camera Link: 100 USB 3.0: 40
Read noise (e ⁻)	0.9 @ 200 MHz 1.2 @ 560 MHz
Sensor cooling 0%	C (up to +30°C ambient)
QE max	60%

Features and Benefits

Compact and light	Highest QE sCMOS sensor
Rolling and true Global (Snapshot) shutter	Maximum exposure and rea freeze frame capture of fast
Industry fastest frame rates	100 fps sustained via Came
FCS Mode	Offers a market-leading spe the temporal demands asso
Low fan vibration	Designed with vibration sen
5.5 megapixel sensor format and 6.5 μm pixels	Delivers extremely sharp res pathology, high content scre
Dual-gain amplifiers	Maximum well depth and lo
12-bit and 16-bit modes	12-bit for smaller file size an
GPU Express	Simplify and optimize data t accelerated GPU processin
Sub-microsecond inter-frame gap	Global Shutter offers down
Water cooled option	Access lowest possible vibr
Extensive FPGA on-head data processing	Essential to ensure best ima
Hardware timestamp	FPGA generated timestamp
Dynamic baseline clamp	Essential to ensure quantita kinetic series.
Spurious noise filter	Real time FPGA filter that id

Key Applications	Zyla 5.5 QE curve
Cell Motility	70
Stem Cells imaging	
Localization microscopy	60
Flow and turbulence studies	
Microfluidics	50
Fluorescence Correlation Microscopy	
Light sheet microscopy	»» × 40
Vesicle Trafficking	cienci
Ion Signaling	ill 30
High throughput screening	Quantum efficiency (%)
	ő



or, providing optimal signal to noise in low light applications.

eadout flexibility across all applications. Snapshot for 'Interline CCD mode' st moving/changing events.

nera Link (full frame). Industry fastest USB 3.0 frame rates.

beed of 27,057 fps from a 2560(h) x 8(v) region of interest, ideally matched to sociated with Fluorescence Correlation Spectroscopy.

nsitive experiments in mind, such as super-resolution microscopy.

esolution over a 22 mm diagonal field of view: Ideal for cell microscopy, digital reening and astronomy.

lowest noise simultaneously, affording extended dynamic range of 33,000:1.

nd absolute fastest frame rates through USB 3.0; 16-bit for full dynamic range.

transfers from camera to Graphical Processing Unit (GPU) card to facilitate ng as part of the acquisition pipeline.

n to 100 ns inter-frame gap, ideal for PIV applications.

pration and -10°C cooling.


nage quality and quantitative fidelity from sCMOS technology.

np with 25 ns accuracy.

300

ative accuracy across the image area and between successive images of a


dentifies and compensates for spurious high noise pixels.



7

Zyla 4.2 PLUS sCMOS

High and broad QE workhorse imaging solution



Zyla 4.2 PLUS is the latest sCMOS technology advancement from Andor. Offering the highest and broadest sCMOS QE profile available, maximizing at 82%, and ideally suited to a wide range of common fluorophores.

The Zyla 4.2 PLUS is also uniquely speedoptimized to deliver a sustained 53 fps at full resolution through a convenient USB 3.0 interface, 77% faster than competing sCMOS cameras. New on-camera intelligence delivers a significant linearity improvement, providing unparalleled quantitative measurement accuracy across the full dynamic range.

The high resolution 4.2 megapixel camera delivers QE performance that is 10% higher than the previous gold-standard sCMOS sensor, coupled with < 1 e- read noise and 33,000:1 dynamic range in a light, compact, low vibration design, intended for both research and OEM usage. The Zyla 4.2 PLUS also includes application specific capability, such as FCS Mode. FCS Mode offers a market-leading speed of 26,041 fps from a 2048(h) x 8(v) region of interest, ideally matched to the temporal demands associated with Fluorescence Correlation Spectroscopy.

New GPU Express facilitates real time data processing.

The Zyla 4.2 PLUS differs fundamentally from Zyla 5.5 in terms of shutter flexibility. Whereas Zyla 5.5 offers both Rolling and true Global shutter modes, the Zyla 4.2 PLUS operates in Rolling shutter mode. However, a mechanism called 'Simulated Global Exposure' is available, whereby a TTL output from the camera can be used to activate a pulsed light source, emulating the Global shutter exposure condition, albeit with less efficiency.

Key Specifications


Active pixels		2048 x 2048
Pixel size (w x h; µm))	6.5 x 6.5
Image area (mm)		13.3 x 13.3
Active area pixel wel	l depth (e)	30,000
Max readout rate (M	Hz)	540
Frame rates (fps)		Camera Link: 100 USB 3.0: 53
Read noise (e ⁻)		0.9 @ 216 MHz 1.1 @ 540 MHz
Sensor cooling	0°C (up t	to +27°C ambient)
QE max		82%

Features and Benefits

82% QE	Highest available photon ca
Market leading USB 3.0 speed	Superb USB 3.0 data trans full resolution, 77% faster th resolution.
FCS Mode	Offers a market-leading spe the temporal demands asso
LightScan PLUS	Reduce background and im to maximise signal and com and Line Scanning Confoca
Low fan vibration	Designed with vibration sen
4.2 megapixel sensor format and 6.5 µm pixels	Delivers extremely sharp res astronomy.
Rolling shutter and simulated Global Exposure	Rolling shutter mode optimi possibility of Rolling shutter
Dual-gain amplifiers	Maximum well depth and lo
> 99.8 % linearity (> 99.9% in low light range)	Unparalleled quantitative ac
GPU Express	Simplify and optimize data t accelerated GPU processin
Dark Noise Suppression (DNS) technology	Extremely competitive low of across range of exposure of
Water cooled option	Access lowest possible vibr
Hardware timestamp	FPGA generated timestamp
Dynamic baseline clamp	Essential to ensure quantita kinetic series.
Spurious noise filter	Real time FPGA filter that id
iCam	Market leading exposure sv

Key Applications Light sheet microscopy Cell Biology Neuroimaging

Neuroimaging
Photoelastic Ptychography
Plant fluorescence
TIRF
FRET
Vesicle Trafficking
Ion Signaling
Microfluidics
C. elegans behavioral mapping
Co-imaging with AFM





apture efficiency across visible/NIR, optimized for all common fluorophores.

sfer efficiency and Zyla's unique 12-bit high speed mode deliver up to 53 fps than competing sCMOS. Follow dynamic processes with improved temporal

beed of 26,041 fps from a 2048(h) \times 8(v) region of interest, ideally matched to sociated with Fluorescence Correlation Spectroscopy.

mprove contrast and resolution in scattering samples. Designed to allow users nfocality concurrently in applications such as Scanned Light Sheet Microscopy and Microscopy.

nsitive experiments in mind, such as super-resolution microscopy.

esolution over a 18.8 mm diagonal field of view; ideal for cell microscopy and

nizes read noise and frame rate. Employ simulated Global shutter method if r spatial distortion.

lowest noise simultaneously, affording extended dynamic range of 33,000:1.

ccuracy of measurement across the full dynamic range

transfers from camera to Graphical Processing Unit (GPU) card to facilitate ng as part of the acquisition pipeline.

dark current of 0.10 $\ensuremath{\text{e}}\xspace$ with fan cooling. Maintains low noise advantage conditions.


pration and -10°C cooling.

np with 25 ns accuracy.

ative accuracy across the image area and between successive images of a


dentifies and compensates for spurious high noise pixels.

witching with minimal overheads.



Zyla - Biologist's Choice

Zyla - Physicist's Choice



Zyla sCMOS has proven a superb camera choice - a price/performance 'workhorse' imager with which to replace or upgrade their existing CCDs including interline cameras.

Example Areas of Application

Lightsheet Microscopy

LightScan PLUS is a new Light Sheet feature set available on Zyla 4.2 PLUS and it has the following benefits:

- Optimize signal to noise AND confocality concurrently
- CycleMax Maximum frame rates with reduced dead-time, no need to reset light sheet for each alternate frame

TIRF Microscopy

The Zyla's fine pixel resolution, great sensitivity, large field of view and fast imaging speed offers a superb choice of platform for following/tracking fast processes at the cell membrane. Multi wavelength TIRF may benefit from Zyla 5.5 in Global shutter.

Super Resolution Microscopy

The low vibration, high QE, low noise and speed capability of Zyla 4.2 PLUS is well suited to single molecule based 'STORM / PALM' approaches. Note, this should be considered distinct from the general needs of single molecule microscopy, which are best served by back-illuminated EMCCD cameras (see Andor iXon EMCCD range)

Physiology / Ion Imaging

The fast frame rate and excellent sensitivity of Zyla is well suited to ion signalling microscopy, Zvla 4.2 PLUS offers sensitivity at speed, but electrophysiology may require the Global shutter exposure mode of Zyla 5.5 to ensure temporal correlation across the whole image.

Cell Motility

Motile cells can be imaged extremely well by the Zyla sCMOS. Its Rolling shutter is suitable for such imaging, but care must be taken of distortive effects if the cell is moving fast. It has been noted that the Global shutter mode was required in Zyla 5.5 to image highly motile sperm cells.

Quality, Throughput and Performance

High Sensitivity and Wide Dynamic Range quantify very weak and very bright structures with one image.


Superb Image Quality - large FOV and high resolution and uniform backgrounds for publication-quality imaging.

Blazingly Fast - studies of cell processes require ever greater temporal resolution.

GPU Express for real time data processing.

Ease of use - designed to get you up and imaging in no time

Flexible - fast or slow, big or small, weak or bright ... Zyla is adaptable all of your imaging



Zyla sCMOS has become a well established detector amongst physicists, biophysicists and astronomers; the advanced combination of speed, sensitivity and dynamic range enabling new ground to be broken.

Example Areas of Application

Lucky / Speckle Imaging

Zyla's fast frame rate and large field of view are ideal for this resolution enhancing technique. GPU Express for real time data processing.

Solar Astronomy

Fast frame rates, wide dynamic range and great linearity present a very formidable solution to the specific detector needs of next generation large solar telescopes.

Adaptive Optics

Accessing > 1000 fps using ROIs renders the Zyla an ideal Wavefront detector. Use with a data splitter to enable direct data access. GPU Express for real time data processing.

Bose Einstein Condensation

The QE profile of Zyla PLUS is very good in the red/NIR region, ideal for BEC of Rb.

Particle Imaging Velocimetry (PIV)

The true Global Shutter mode of Zyla 5.5 facilitates an inter-frame gap of down to 100 ns.

Fluorescence Correlation Spectroscopy The FCS mode, now available on Zyla 4.2 PLUS provides market leading speeds of 26,041 fps from a 2048(h)x8(v) region of interest.

X-Ray / Neutron Tomography The Zyla can be readily lens coupled to scintillators and phosphors, presenting a high resolution, sensitive and fast solution for tomography.*



coupled camera, superb for fast indirect X-Ray applications such as tomography or non-destructive testing.



Fluorescence Correlation Spectroscopy



Check out Andor's new Zyla HF Fiber Optic

Performance and Adaptability

Dual Amplifier - access lowest read noise and full well depth simultaneously.

1,000 fps - extremely fast frame rates through user definable Region of Interest control. GPU Express for real time data processing.

Global shutter - this important mode completely avoids spatial distortion, and ensures temporal correlation across all regions of the sensor. Sub-microsecond interframe gap - ideal for PIV applications.

Low dark current - extremely low dark current ensuring minimized hot pixel blemishes.

Cooling options - standard camera air cools to 0°C up to +35°C ambient. Water cooled option is available

Blemish correction maps and advanced control are available upon request

Compact and Light - the extremely small volume footprint of Zyla renders it adaptable to intricate optical set-ups.

Upgrade your performance

Zyla sCMOS - Disruptive technology, familiar price.

What was the benchmark imaging detector?

Interline CCD technology has been the dominant 'workhorse' detector type for fluorescence cell microscopy for almost 15 years, the benchmark sensor being a 1.4 megapixel, 6.45 µm pixel size device, offering typical read noise floor between 5 and 8 e- rms at modest frame rates of 11-12 fps.

Upgrade your microscope using the imaging superiority of Zyla sCMOS:

Fundamentally, sCMOS technology has been conceptualized as a vastly superior alternative to Interline CCDs. Indeed, Andor's Zyla sCMOS offers dramatically higher performance, yet remains within the same price bracket as Interline cameras, and is ideally placed to become the new gold standard 'workhorse' laboratory detector. Importantly, Zyla uniquely comes with both Rolling and true Global (Snapshot) shutter.

In particular, Global shutter offers a simple Snapshot imaging capability, directly analogous to that of Interline CCDs, offering zero image distortion and perfect for synchronizing to peripheral devices.

Benefits of upgrading...





10x 16x more sensitive




The table below compares some typical performance specifications of a 1.4MP Interline CCD to those of the new Zyla sCMOS, outlining the approximate factor improvement that is available.

Parameter	Typical Interline CCD	Zyla 5.5 sCMOS	sCMOS Factor Improvement
Read noise	6 e-	0.9 e-	5x more sensitive
Sustained frame rate	12 fps @ 1.4 MP	100 fps @ 5.5 Megapixel* 200 fps @ 1.4 Megapixel ROI*	16x faster
Dynamic range	2,250:1	33,000:1	10x more dynamic range
Sensor format	1.4 Megapixel	5.5 Megapixel	4x more pixels

* Frame rates provided are for the '10-tap' Cameralink model. An even more affordable USB 3.0 model is available, in the case of the Zyla 5.5 offering 40 fps sustained at full 5.5 megapixel resolution.

Why is Zyla such an impressive all-round imager?

- Superior performance vastly superior to Interline across key performance parameters. • Rolling and Global exposures (Zyla 5.5) - Zyla is unique in offering both these exposure modes in one camera. Global shutter (Snapshot) is directly analogous to the Interline
- exposure mechanism. • Image quality - a huge amount of effort and on-camera (FPGA) intelligence has gone into optimizing image quality in Zyla.
- Flexibility fast, slow, weak, bright, pixel binned, region of interest, Rolling shutter, Snapshot shutter... Zyla is adaptable to a broad gamut of application requirements.
- Affordable price with so many superb features, we have endeavoured to make the Zyla accessible to every lab. Request a quote, you'll be pleasantly surprised!
- USB 3.0 'Plug and play' interface with industry fastest frame rates.



Key Biological Applications

Rey biological Applications	
Live cell imaging	
Widefield fluorescence microscopy	
Developmental biology	
Embryo studies	
Physiology / Ion Imaging	
Neuroscience / Vesicle Transport	
Parasitology	
Cell Motility	





Extended Dynamic Range

The Andor Neo and Zyla cameras are designed to make use of the innovative dual 'column-level' amplifier design of the sensors.

Traditionally, sensors require that the user must select up-front between high or low amplifier gain (i.e. sensitivity) settings, depending on whether they want to optimize for low noise or maximum well depth. The dual amplifier architecture of the sCMOS sensor circumvents this need, in that signal can be sampled simultaneously by both high and low gain amplifiers. As such, the lowest noise of the chip can be harnessed alongside the maximum well depth, affording widest possible dynamic range of up to 33,000:1.



Zyla uniquely offers both 12-bit and 16-bit modes. **12-bit** for smaller file size & absolute fastest frame rates through USB 3.0; 16-bit for full dynamic range.

Dual Amplifier Architecture


Each column within each half of the sensor is equipped with dual column level amplifiers and dual analog-to-digital converters (ADC).

This architecture was designed to simultaneously minimize read noise and maximize dynamic range. The dual column level amplifier/ADC pairs have independent gain settings, and the final image is reconstructed by combining pixel readings from both the high gain and low gain readout channels to achieve an unprecedented intra-scene dynamic range from the relatively small 6.5 µm pixel pitch.

High contrast image recorded with dual amplifier 16-bit mode of Neo

Pixels sampled by high gain amplifier (~800 counts)

Pixels sampled by low gain amplifier (~8,000 counts)





the 'Dual Amplifier Dynamic Range' technical note

Lowest Noise Floor

Andor's ultra sensitive sCMOS cameras have broken new ground in offering down to 0.9 electron rms read noise, without signal amplification technology.

What is truly extraordinary is that this performance level is achievable at 30 fps, representing 200 MHz pixel readout speed. Furthermore, even at full readout speed, the read noise floor is negligibly compromised, maintaining down to 1.3 e- rms at 100 fps. For the best CCD cameras to even approach 2 electrons noise, a readout speed of 1 MHz or slower is required. This minimal detection limit renders Andor's sCMOS cameras suitable for a wide variety of challenging low light imaging applications.

Spurious Noise Filter

Neo and Zyla platforms both come equipped with an optional in-built FPGA filter that operates in realtime to reduce the frequency of the occurrence of high noise pixels that would otherwise appear as spurious 'salt and pepper' noise spikes in the image background.




Neo sCMOS

Comparative low light images taken with Neo sCMOS (1.3 electrons read noise @ 560 MHz) vs. Interline CCD (5 electrons read noise @ 20 MHz), displayed with same relative intensity scaling.

technical note



Readout Speed (MHz)	Neo 5.5 Readout No	ise (e [.])
	Rolling Shutter	Global Shutter
200	1	2.3
560	1.3	2.5
Readout Speed (MHz)	Zyla 4.2 PLUS Reade	out Noise (e ⁻)
	Rolling Shutter	
200	0.9	
540	1.1	



(a) LED signal in a light-tight imaging enclosure, intensity ~30 photons/pixel; (b) Fluorescently labelled fixed cell using a CSU-X spinning disk confocal microscope (x60 oil objective), each 100 ms exposure, same laser power,

Rapid Frame Rates

The parallel readout nature of sCMOS means it is capable of reaching very rapid frame rates of up to 100 full frames per second, and much faster with region of interest.

Distinctively, this is accomplished without significantly sacrificing read noise performance, markedly distinguishing the technology from CCDs. Andor's sCMOS cameras are uniquely designed to harness this speed potential.

Zyla's speed optimized USB 3.0 interface delivers an unparalleled 40 fps from 16-bit mode and 53 fps from Zyla's unique 12-bit mode (4.2 megapixel array).



Array Size	Zyla 5.5 USB 3.0 Rolling Shutter	Global Shutter	Zyla 5.5 10-tap Rolling Shutter	Global Shutter
2560 x 2160	40(30)	40 (30)	100 (75)	49 (49)
2048 x 2048	53 (40)	52 (39)	105 (98)	52 (52)
1920 x 1080	107 (80)	98 (80)	200 (200)	97 (97)
512 x 512	422 (422)	201 (201)	422 (422)	201 (201)
128 x 128	1,691 (1,691)	716 (716)	1,691 (1,691)	716 (716)
2048 x 8 (FCS mode)	13,020 (10,250)	4,008 (4,008)	27,057 (27,057)	4,008 (4,008)

Array Size	Zyla 4.2 PLUS 10-tap Rolling Shutter	Zyla 4.2 PLUS USB 3.0 Rolling Shutter
2560 x 2160	-	-
2048 x 2048	101 (101)	53 (40)
1920 x 1080	192 (192)	107 (80)
512 x 512	406 (406)	406 (406)
128 x 128	1,627 (1,627)	1,627 (1,627)
2048 x 8 (FCS mode)	26,041 (26,041)	13,020 (10,250)

Maximum frame rates achievable from the Zyla 5.5 and Zyla 4.2 PLUS sCMOS USB 3.0 and 10-tap Camera Link versions 12-bit (16-bit)

iCam Fast Exposure Switching

Neo and Zyla benefit from Andor's iCam technology, an innovation that ensures minimal overheads associated with fast exposure switching.

This is particularly important during multi-color microscopy acquisition protocols, whereby it is necessary to repeatedly and rapidly flip between pre-set exposure times matched to the relative signal intensity of each fluorophore.

iCam offers market leading acquisition efficiency, whether software or externally triggered.


GPU Express

The Andor GPU Express library has been created to simplify and optimize data transfers from camera to a CUDAenabled NVidia Graphical Processing Unit (GPU) card to facilitate accelerated GPU processing as part of the acquisition pipeline. Integrates easily with Andor SDK3.

Data Flow Monitor

The sCMOS sensor in Neo and Zyla is capable of extremely fast data read rates, but this in itself imposes considerable challenges.

For sustained kinetic series measurements it is possible to be rate limited by: (a) bandwidth of the Camera Link interface connecting the camera to the PC (b) hard drive write speed



The Data Flow Monitor, available through Andor Solis acquisition and analysis software, has been innovated to provide a simple visual tool that enables you to instantly ascertain if your acquisition parameters will result in a rate of data transfer that is too fast for either interface or hard drive. It will also determine if the kinetic series size is within the capacity of camera memory, hard drive space or PC RAM.




of Camera Link data transfer bandwidth and Hard Disk Drive write speed.



e.g. 2 - Hard Disk Drive will not write data fast enough for the requested kinetic series. Advised to first reduce data rate.





e.g. 1 - Requested kinetic series within capability



Our experiments with Andor's new sCMOS camera have been highly encouraging. The combination of very low noise sensitivity at rapid frame rates, coupled with high pixel resolution and large dynamic range, will enable us to investigate single molecules at timescales which were previously not accessible.

Prof Stefan Diez, Heisenberg Professorship for BioNanoTools, Max Planck Institute of Molecular Cell Biology and Genetics, Dresden, Germany

Rolling and Global (Snapshot) Shutter Modes

Neo 5.5 and Zyla 5.5 are distinct in offering both Rolling shutter and true Global shutter modes from the same sensor, such that the most appropriate mode can be selected dependent on application requirements.

The 5.5 megapixel sCMOS sensor that is in Zyla 5.5 and Neo 5.5 sCMOS cameras uniquely offers both Rolling and Global exposure modes. This provides superior application and synchronization flexibility and the ability, through Global exposure, to closely emulate the familiar 'Snapshot' exposure mechanism of Interline CCDs.


Rolling and Global Shutter Mechanisms

Rolling and true Global shutter modes describe two distinct types of exposure and readout sequence.


In Rolling shutter, available in all Andor sCMOS cameras, different lines of the array are exposed at different times as the read out 'wave' sweeps through the sensor. 10 ms are required at the start to 'activate' the sensor to expose, and then 10 ms are required at the end to readout the sensor. Use when not synchronizing to peripheral devices and only when there is a minimal risk of spatial distortion from slow moving sample.

In true Global shutter, offered in both Neo 5.5 and Zyla 5.5 models, each pixel in the sensor begins the exposure simultaneously and ends the exposure simultaneously. This provides a true 'Snapshot' exposure capability for moving samples that is both 'photon-efficient' and easy to synchronize to, especially useful for 3D / 4D microscopy. Zyla 4.2 PLUS, while utilizing a Rolling

Global shutter exposure and readout (single scan)



Rolling shutter exposure and readout (single scan)



shutter sensor, offers a Simulated Global Exposure mechanism to overcome risk of spatial distortion. This mechanism is more elaborate and less photon/time efficient than true Global shutter. True Global Shutter is also essential for applications such as Particle Imaging Velocimetry (PIV), where sub-microsecond inter-frame gaps are required.



What should I be aware of as a buyer? Beware of 'Gen II' claims!

This topic carries particular relevance, as not all 'scientific CMOS' cameras on the market offer a choice of Rolling and true Global exposure. Most offer one or the other. In fact, a sensor that is currently being widely positioned in the market as "Gen II", is actually the sensor used in the Zyla 4.2 PLUS. However, we would not go so far as to describe this as a 2nd generation sensor, as it achieves a higher Quantum Efficiency at the notable expense of true Global shutter capability.

The customer is highly advised to make an informed choice before purchasing, since in microscopy the more photon-efficient Global shutter approach can actually result in a higher signal to noise, faster synchronized frame rates and non-distorted images.

Global Exposure is Distortion Free

If light is falling on the sensor during the 'transient' phases (first 10 ms and final 10 ms) of the Rolling shutter exposure mechanism, and an object is moving during this time, then there is a chance of some degree of spatial distortion. The degree of distortion is dictated by the relative size, direction and speed of the object. Global shutter avoids spatial distortion since there are no 'transient' exposure phases.

These images show the head of a moving sperm cell, imaged by the Neo sCMOS in both Global and Rolling shutter. Distortion of the shape of the sperm head is evident in Rolling shutter.

Key Benefits of true Global shutter (Zyla 5.5 and Neo 5.5)

Global shutter in particular is viewed as an important mode for the biologist, as its benefits are deeply synergistic with the core imaging requirements of live cell microscopy.

No spatial distortion - avoiding the spatial distortion risk of Rolling exposure

3D / 4D microscopy - recommended for synching to peripheral switching devices

Higher signal to noise due to reduced dead time offering higher 'effective' QE

Simplicity - all the benefits of an 'Interline exposure mode'

Compatible with continuous or pulsed light sources

Sub-microsecond inter-frame gap, ideal for PIV applications






To find out more about how to synchronize to Rolling and Global shutter (and to view our FAQs) visit andor.com/learning and read the 'Synchronizing to Andor sCMOS Cameras' technical note.



Rolling



Correct

Distorted

For our work on quantifying red blood cell velocity in the retinal capillaries we elected to operate in Global shutter mode, which produces minimally distorted images. When operating in Rolling shutter mode we observed significant image warping, even for moderate eve movements.

Dr. Phillip Bedggood, Metha Laboratory, Department of Optometry and Vision Sciences, University of Melbourne, Australia


Deep Thermoelectric Cooling


Andor's Neo offers the deepest sensor cooling available from any CMOS imaging camera on the market, minimizing both dark current and hot pixel blemishes. Additionally, through the use of water cooling, the fan can be switched off in the software to minimize camera vibration; ideal for set-ups that are particularly vibration sensitive.

Neo Cooling Temperature	Dark current (e ⁻ /pix/sec)
-30°C (fan cooling)	0.015
-40°C (10°C liquid)	0.007



The temperature sensor in the **Neo** and **Zyla sCMOS** cameras measures with a **thermostatic precision of 0.05°C.**





Thermal noise can sacrifice the sCMOS low detection limit. Low light images recorded with a Neo sCMOS camera at +5°C and -40°C sensor cooling temperatures; 50 sec exposure time; 200 MHz readout giving 1 electron read noise.

Deep TE cooling is useful for a number of reasons:


Minimization of dark current

sCMOS cannot be considered a truly flexible, workhorse camera unless dark current contribution has been minimized. Deep cooling means the low noise advantage can be maintained under all exposure conditions.



Minimization of hot pixel blemishes

Hot pixels are spurious pixels with significantly higher dark current than the average and can be problematic even under relatively short exposure times. Cooling has a major influence in minimizing the occurrence of such events, offering both an aesthetically cleaner image and a greater number of unfiltered, usable pixels.



Hot pixel blemishes are significantly reduced at deeper cooling temperatures, requiring much reduced pixel correction. Uncorrected images are shown above for 1 sec exposure.

Minimization of vibration

Many optical configurations are sensitive to vibrations from the camera fan. Andor's Neo offers:

(a) Two fan speeds
(b) The ability to turn off the fan, either temporarily or permanently if flowing liquid through the camera (the latter also allows the Neo to be stabilized at -40°C)

UltraVac[™] (Neo only)

The Andor Neo is the only vacuum housed CMOS sensor available on the market, offering superior quality, performance and longevity.

Andor's proprietary UltraVac[™] process has a proven track record of field reliability, accumulated over more than 15 years of shipping high-end vacuum cameras. Using a proprietary technique, we have adapted these process for use with the additional connections associated with the sCMOS sensor.

Key Features

Permanent hermetic vacuum seal
Sustained deep TE cooling
No maintenance / re-pumping
No risk of condensation
Minimize out-gassing




5 Year Vacuum Warranty

Our faith in the unique sCMOS vacuum process used in Neo means that we are proud to offer an extensive 5 year warranty on the vacuum enclosure.







Schematic of the Neo sCMOS permanent vacuum head





Advanced FPGA On-head Processing

Andor's Neo and Zyla cameras are each equipped with considerable FPGA processing power. This is essential in order to dynamically normalize data at the pixel level for minor variations in bias offset, thus eradicating fixed pattern noise associated with this CMOS phenomenon.

This superior dynamic processing capability is also utilized to optionally filter the small percentage of spurious noise pixels from the image.

CMOS data requires compensation for fixed pattern variation. This is accomplished in real time for every pixel within the FPGA of Andor's sCMOS cameras, essentially eliminating this noise source from the image.

Pixel-level bias offset compensation

The advanced processing power and memory capacity permits implementation of bias offset compensation for every pixel in the array. This ultimately relates to a lower noise background.

Dynamic baseline clamp

Key Features

and well depth

noise

loss

Superb 30,000 electron well depth

With pixel compensation

A real time algorithm that uses dark reference pixels on each row to stabilize the baseline (bias) offset. Necessary to ensure quantitative accuracy across each image and between successive images.

Spurious noise filter


Without pixel compensation

An optional real time filter that identifies and compensates for 'spurious' high noise pixels that are greater than 5 electrons (< 1% of all pixels)

Andor offer the capability to switch off interpolative pixel filtering and provision of custom blemish maps, important for applications such as super-resolution microscopy or astronomy.

6.5 µm pixel size combined with 30,000 electron well depth

The 6.5 µm pixel present in Neo and Zyla has been specifically designed to offer an optimal balance of optical resolution, photon collection area and well depth. This pixel size has been determined to provide ideal over-sampling of the diffraction limit in typical cell microscopy with x 60 and x 100 objectives.



Dynamic Baseline Clamp affords superior background image quality (Dark image, Rolling shutter, 10 ms exposure)

Large Field of View

The multi-megapixel sensors present in the Neo and Zyla offer an extended field of view, markedly exceeding the FOV available from alternative Interline CCD devices.

Flexibility is key however, and if a large FOV is not required for a particular application, Neo and Zyla offer a range of pre-selected ROI sizes at the click of a button, as well as user defined (with single pixel granularity).

Key Features

21.8 mm diagonal (Zyla 5.5 and Neo 5.5); 18.8 mm diagonal (Zyla 4.2 PLUS) Closely matched to modern microscopes Pre-selected ROIs to quickly opt for smaller FOV if required x 3.5 larger than popular 512 x 512 EMCCD sensor x 3.9 larger than popular 1.4 MP Interline CCD

sensor

Combine large FOV with dual wavelength image splitter (OptoSplit)

Field of View Comparison Neo 5.5 and Zyla 5.5 Field of View vs popular 1.4 megapixel Interline CCD

Trigger Mode

Time Lapse

Continuous

External Exposure

External Start


Comprehensive trigger functionality


Neo and Zyla offer a selection of advanced trigger modes, designed to provide tight synchronization of the camera within a variety of experimental set-ups. Triggering is compatible with both Rolling and Global shutter modes.

• External TTL, Software and Internal trigger (including Simulated Global Exposure - Zyla 4.2 PLUS)


• 'Time Lapse' and 'Continuous' (overlapped) kinetic series

• Fast exposure switching (iCam)











Description

- Each exposure started by a trigger event (e.g. TTL rising edge). Exposure duration is internally defined.
- Exposures run back to back with no time delay between them. Exposure time defined by time between consecutive trigger events.
- Exposure time defined by TTL width (sometimes known as 'bulb mode').
- TTL rising edge triggers start of internally defined kinetic series.
- Available Neo and Zyla trigger modes, applicable to both Rolling and Global shutter.

Hyperspectral Imaging of Fluorescent Beads. Image courtesy of Dr. Ethan Schonbrun.

Software Solutions

Andor Solis

Solis is a ready to run Windows package with rich functionality for data acquisition and image analysis/processing. Available on 32-bit and 64-bit versions of Windows (7, 8 and 10). Andor Basic provides macro language control of data acquisition, processing, display and export.

Andor iQ

A comprehensive multi-dimensional imaging software package. Offers tight synchronization of EMCCD with a comprehensive range of microscopy hardware, along with comprehensive rendering and analysis functionality. Modular architecture for best price/performance package on the market.



Andor SDK

A software development kit that allows you to control the Andor range of cameras from your own application. Available as 32 and 64-bit libraries for Windows (7, 8 and 10) and Linux. Compatible with C/C++, C#, Delphi, VB6, VB.NET, LabView and Matlab.

Bitplane Imaris®

Imaris delivers all the necessary functionality for visualization, segmentation and interpretation of multidimensional datasets. By combining speed, precision and intuitive ease-of-use, Imaris provides a complete set of features for handling multi-channel image sets of any size up to 50 gigabytes.

Third party software compatibility
Drivers are available so that the iXon range

can be operated through a large variety of third party imaging packages, including:

Metamorph (Molecular Devices Corporation)
NIS Elements (Nikon)
LAS (Leica)
Xcellence (Olympus)
Image Pro (Media Cybernetics)
MicroManager (UCSF)
Till Photonics Live Acquisition (Till Photonics)
Imaging Workbench (Indec)
WinFluor (University of Strathclyde)
Maxim DL (Diffraction Limited)

Maxim DL (Diffraction Limited) LabView (National Instruments)

Matlab (MathWorks)

Extensive Imaging Portfolio

The Andor Imaging Range

Have you found what you are looking for? As an alternative to the sCMOS cameras, Andor offers an extensive portfolio of high performance low light imaging camera technologies.



iXon EMCCD High performance EMCCD platform

Single photon sensitive and back
Industry fastest frame rates
-100°C cooling
Flexible yet intuitive
Quantify in electrons or photons



Zyla sCMOS Fast, sensitive, compact, light sCMOS < 1 electron read noise @ 30 fps 5.5 and 4.2 Megapixel sensors / 6.5 µm QE up to 82% 0°C cooling at +27°C ambient 100 fps sustained (10-tap Camera Link) Cost effective USB 3.0 option 16-bit data range



Neo sCMOS Vacuum cooled, lowest noise sCMOS 1 electron read noise @ 30 fps 5.5 Megapixel / 6.5 µm -40°C vacuum cooling 30 fps sustained; 100 fps burst

4 GB on head memory 16-bit data range Fan off vibration free mode





and back-illuminated

25

Notes

Application & Technical Notes

The following section is dedicated to providing a greater depth of understanding of the performance and innovations associated with the Andor scientific CMOS camera platform. Deeper insight is provided into areas such as the unique dual amplifier architecture (for extended dynamic range), sCMOS read noise distribution, dark noise effects, vacuum sensor protection and Rolling vs. Global shutter readout modes.

We also present a comprehensive overview of how new sCMOS technology compares to existing 'gold standard' scientific imaging cameras such as Interline CCD and EMCCD technology.

Visit the Learning Center now to discover more at andor.com/learning.

Key sCMOS resources include:

- LightScan PLUS
- Rolling and Global Shutter
- Dual Amplifier Dynamic Range
- The Importance of TE Cooling to sCMOS Technology
- Comparing sCMOS With Other Scientific Detectors
- Andor sCMOS PC Recommendations and Data Flow Considerations
- Understanding Read Noise in sCMOS
- Interpolative Blemish Corrections on sCMOS
- PIV Mode for Neo and Zyla







Customer Support

Andor products are regularly used in critical applications and we can provide a variety of customer support services to maximize the return on your investment and ensure that your product continues to operate at its optimum performance.

Andor has customer support teams located across North America, Asia and Europe, allowing us to provide local technical assistance and advice. Requests for support can be made at any time by contacting our technical support team at andor.com/support.

Andor offers a variety of support under the following format:

- On-site product specialists can assist you with the installation and commissioning of your chosen product
- Training services can be provided on-site or remotely via the Internet
- A testing service to confirm the integrity and optimize the performance of existing equipment in the field is also available on request.

A range of extended warranty packages are available for Andor products giving you the flexibility to choose one appropriate for your needs. These warranties allow you to obtain additional levels of service and include both on-site and remote support options, and may be purchased on a multi-year basis allowing users to fix their support costs over the operating life cycle of the products.



Head Office

7 Millennium Way Springvale Business Park Belfast BT12 7AL Northern Ireland Tel: +44 (0)28 9023 7126 Fax: +44 (0)28 9031 0792

North America

300 Baker Avenue Suite 150 Concord, MA 01742 USA Tel: +1 860-290-9211 Fax: +1 860-290-9566

Japan

5F IS Building 3-32-42 Higashi-Shinagawa Tokyo 140-0002 Japan Tel: +81-(0)3-6732-8968 Fax: +81-(0)3-6732-8939

China

Unit 1, Building A, No. 66 Zhufang Road, Haidian District, Beijing 100085 P. R. China Tel: +86 (0)10-8271-9066 Fax: +86 (0)10-8271-9055

Find us on



Andor[®], the Andor logo and iXon[®] are trademarks of Andor Technology Ltd. Changes are periodically made to the product and specifications are subject to change without notice.

