


# iXon Ultra & Life 897

Version 2.3 rev 28 Jan 2025



## Hardware Guide

Covering the iXon Ultra 897 and iXon Life 897 models

© Andor Technology Ltd. 2025

<b>Section 1: Introduction .....</b>	<b>10</b>
1.1 Technical Support .....	11
1.2 Disclaimer.....	12
1.3 Copyright and Protective Notices.....	12
1.4 Trademarks and Patent Information .....	12
1.5 Supplied Components.....	13
1.5.1 Optional Components .....	13
1.6 Camera Power Supply Unit (PSU) .....	14
1.6.1 Working with Electronics.....	14
1.7 Prevention of Condensation .....	15
1.8 EM Gain Ageing.....	15
1.8.1 Minimizing EM Gain Ageing in your iXon Ultra/Life Camera .....	15
1.9 Minimizing Particulate Contamination.....	16
1.10 Software .....	16
<b>Section 2: Product Overview .....</b>	<b>17</b>
2.1 iXon Ultra and Life 897.....	17
2.2 Power and Signal Connections.....	18
2.3 Camera Link (iXon Ultra only).....	19
2.4 External I/O.....	19
2.4.1 Additional Cables (iXon Ultra Only).....	20
2.5 Signal Diagrams .....	21
2.5.1 iXon Ultra and Life Input & Output Timing Hardware.....	21
2.5.2 External Trigger Input (at connector) .....	21

<b>Section 3: Installing the iXon Life and Ultra 897 .....</b>	<b>22</b>
3.1 PC Requirements .....	22
3.2 Connecting the Camera.....	22
3.3 Cooling the CCD .....	22
3.3.1 Sources of Heat Generation .....	23
3.4 Connecting a Cooling System .....	23
3.4.1 Hose Connections.....	23
3.4.2 Coolant Recommendations.....	23
3.4.3 Connecting the Coolant Hoses .....	24
3.4.4 Removing the Coolant Hoses.....	24
3.5 Mounting Posts.....	25
3.5.1 Attaching Mounting Posts to the Camera .....	25
3.6 Installing Andor Solis Software for Windows 10 .....	26
3.7 New Hardware Wizard .....	26
3.8 Start-up Dialog.....	27
<b>Section 4: Features and Functionality .....</b>	<b>28</b>
4.1 EMCCD Operation .....	28
4.1.1 Structure of an EMCCD.....	28
4.1.2 EM Gain & Read Noise .....	30
4.1.3 EM Gain ON vs EM Gain OFF .....	31
4.1.4 Multiplicative Noise Factor and Photon Counting.....	32
4.1.5 EM Gain Dependence and Stability .....	33
4.1.6 RealGain™: Real and Linear gain .....	34
4.1.7 EM Gain Ageing: What causes it and how is it countered?.....	35
4.1.8 Gain and Signal Restrictions .....	35
4.1.9 EMCAL™.....	36
4.2 Sensor Cooling Considerations .....	37
4.2.1 Thermoelectric Cooling.....	37
4.2.2 Readout Mode and Heat Generation in the EMCCD .....	37
4.2.3 Heatsink "hot side" Temperature.....	37
4.2.4 Cooling Options .....	38
4.2.5 Adjusting the Sensor Operating Temperature.....	38
4.2.6 Cooling Fan Settings.....	39

<b>4.3</b>	<b>Sensor Readout Optimisation.....</b>	<b>39</b>
4.3.1	Sensor Pre-amp Options.....	39
4.3.2	Variable Horizontal Readout Rate.....	40
4.3.3	Variable Vertical Shift Speed.....	41
4.3.4	Output Amplifier Selection .....	42
4.3.5	Baseline Optimization .....	42
4.3.5.1	Baseline Clamp.....	42
4.3.6	Binning and Sub Image options.....	43
<b>4.4</b>	<b>Acquisition Options.....</b>	<b>45</b>
4.4.1	Capture Sequence in Frame Transfer (FT) Mode.....	45
4.4.1.1	Points to consider when using FT Mode.....	46
4.4.2	Capture Sequence in Non-Frame Transfer Mode (NFT) with an FT EMCCD.....	47
4.4.2.1	Points to note about using an FT EMCCD as a standard EMCCD.....	48
4.4.3	Capture Sequence for Fast Kinetics (FK) with an FT EMCCD .....	48
4.4.3.1	Points to consider when using Fast Kinetics mode.....	49
	<b>Section 5: Triggering Information .....</b>	<b>50</b>
<b>5.1</b>	<b>Keep Clean Cycles.....</b>	<b>50</b>
5.1.1	Idle Keep Clean Cycle.....	50
5.1.2	Internal Keep Clean Cycle .....	51
5.1.3	External Keep Clean Cycle.....	52
<b>5.2</b>	<b>Triggering Modes.....</b>	<b>53</b>
<b>5.3</b>	<b>Triggering Options in Frame Transfer (FT) Mode .....</b>	<b>54</b>
5.3.1	Internal Triggering.....	54
5.3.2	External Triggering.....	55
5.3.3	iCam Enhanced Trigger Mode .....	56
<b>5.4</b>	<b>Triggering Options in Non-Frame Transfer (NFT) Mode.....</b>	<b>57</b>
5.4.1	Internal Triggering (NFT).....	57
5.4.2	External & Fast External (NFT) Triggering .....	58
5.4.3	External Exposure (NFT) Triggering .....	59
5.4.4	Software (NFT) Triggering .....	60
<b>5.5</b>	<b>Triggering Options in Fast Kinetics (FK) Mode .....</b>	<b>61</b>
5.5.1	Internal (FK).....	61
5.5.2	External (FK) .....	62
5.5.3	External Start (FK) .....	63
<b>5.6</b>	<b>Shuttering (iXon Ultra Only).....</b>	<b>64</b>

---

<b>Section 6: Additional Functionality</b>	<b>65</b>
6.1 Count Convert	65
6.2 OptAcquire	66
6.2.1 OptAcquire modes	67
6.3 Pushing Frame Rates with Cropped Sensor Modes	68
6.3.1 Cropped Sensor Mode Frame Rates	69
6.4 Advanced Photon Counting in EMCCDs (iXon Ultra Models Only)	70
6.4.1 Photon Counting by Post-Processing	72
6.5 Spurious Noise Filter	73
<b>Section 7: Troubleshooting</b>	<b>74</b>
7.1 Unit does not switch on	74
7.2 Support Device not recognised when plugged into PC	74
7.3 Temperature Trip Alarm sounds (continuous tone)	74
7.4 Camera High FIFO Fill Alarm	75
7.5 Other Possible Causes for Stoppages in Acquisition	75
7.6 Use of Multiple High Speed USB 2.0 I/O on One Camera	76
7.6.1 Mapping each port	76
<b>Section 8: Maintenance</b>	<b>77</b>
8.1 Cleaning Camera Exterior	77
8.2 Regular Checks	77
8.2.1 Annual Electrical Safety Checks	77
8.3 Replacement Parts	77
8.4 Fuse Replacement	77

---

Appendix A: Technical Specifications.....	78
Appendix B: Mechanical Drawings.....	82
Appendix C: Dew Point Graph .....	83
Appendix D: EMCCD Technology.....	84
Appendix E: Glossary.....	88
Appendix F: Other Information.....	94
Appendix G: iXon China RoHS Hazardous Substances Declaration.....	95

## Revision History

Version	Released	Description
1.0	01 Jul 2012	Initial Release
1.1	09 Jan 2013	Version ID added to footer of each page for clearer identification of version and revision Minor updates to formatting and text throughout to improve clarity Section 2.5 "Water Pipe Connectors" updated Section 2.7 "Cooling" updated to provide improved information on cooling Sections 3.2 to 3.2.4 updated to provide improved information on cooling Section 4.8 "Cameralink" updated to clarify Cameralink connection functionality References to obsolete ACZ-03452 updated to ACZ-03463 References to CCD updated where relevant to EMCCD
1.2	17 Sep 2014	Branding updated to latest format throughout. Manual Structure and format updated throughout. Updated Specifications to show iXon is a Class A device Count Convert: Added base mean level for EX and EXF sensor variants to 500 (Section 6.1) Added frame rates for Optically Centred Crop Mode (Section 6.3.1)
1.3	13 May 2015	Further detail provided for power supply requirements (Section 1.6) Removed external exposure for frame transfer mode- not supported (Section 5.3)
1.4	04 Jul 2017	Updated to include iXon Life models.
1.5	29 Nov 2017	Minor edits to align weight, storage temp and clearance specifications across documentation.
1.6	25 Jan 2018	Mechanical drawings updated - water tap dimensions (Appendix A). Base mean level information updated (Section 6.1)
1.7	08 Jan 2018	Frame rate figures updated - matching tests on latest software (Section 6.3.1) Added section 2.4.1 Additional Cables (iXon Ultra Only).
1.8	11 Apr 2019	Updated US and Japan addresses, fixed hyperlinks.
1.9	30 Sep 2019	Updated mechanical drawing of backplate, updated China office address and phone number.
2.0	21 Jun 2022	Removed links to MyAndor, added links to Downloads area. Updated styling throughout. Added note on DAC output availability. Added China RoHS table to appendix. Added Reach statement.
2.1	26 Sep 2023	Replaced Fig 28. Removed support for Windows 8.
2.2	01 Oct 2024	Removed BVF and #EX options due to obsolescence (Appendix A).
2.3	28 Jan 2025	Updated "PC Requirements" to include Windows 11. (Section 3.1)

## Updates to the Manual

Changes are periodically made to the product, and these will be incorporated into new editions of the manual. Please check for new releases of the manual at: [andor.oxinst.com/downloads](https://andor.oxinst.com/downloads). If you find an issue in this manual, please contact your customer support representative (Section 1.1) with a description of the issue.

---

## Safety and Warning Information



### PLEASE READ THIS INFORMATION FIRST

1. To ensure correct and safe operation of this product, please read this guide before use and keep it in a safe place for future reference.
2. If the equipment is used in a manner not specified by Andor, the protection provided by the equipment may be impaired.
3. Before using the system, please follow and adhere to all warnings, safety, manual handling and operating instructions located either on the product or in this guide.
4. This camera is a precision scientific instrument containing fragile components. Always handle with care.
5. Do not expose the product to extreme hot or cold temperatures.
6. Ensure that the ventilation slots in the camera case are free from blockages.
7. Do not expose the product to open flames.
8. Do not allow objects to fall on the product.
9. Do not expose the product to moisture, wet, or spill liquids on the product. Do not store or place liquids on the product. If spillage occurs on the product, switch off power immediately and wipe off with dry, lint-free cloth. If any ingress has occurred or is suspected, unplug mains cable, do not use, and contact Andor service.
10. The product contains components that are extremely sensitive to static electricity and radiated electromagnetic fields, and therefore should not be used, or stored, close to EMI/RFI generators, electrostatic field generators, electromagnetic or radioactive devices, or other similar sources of high energy fields.
11. Electromagnetic Compatibility: This product has been designed and tested to perform successfully in a normal (basic) electromagnetic environment, e.g. a typical life science test laboratory, as per the EU EMC Directive. It is not designed to operate in a harsh electromagnetic environment, e.g. close to the following equipment: EMI/RFI generators, electrostatic field generators, electromagnetic or radioactive devices, plasma sources, arc welders, x-ray instruments, intense pulsed sources, or other similar sources of high energy fields whose emissions are not within the normal range expected under the EU EMC Directive.
12. Electromagnetic Compatibility: This is a Class A product. In a domestic environment this product may cause electromagnetic interference, in which case the user may be required to take adequate measures.
13. This product is not designed to provide protection from ionising radiation. Any customer using this product in such an application should provide their own protection.
14. This product is for use in research laboratories and other controlled scientific environments.
15. This product has not been designed and manufactured for the medical diagnosis of patients.
16. Use only the power supply cord provided with the system for this unit. Should this not be correct for your geographical area, contact your local Andor representative.
17. The camera should be mounted so that the mains supply can be easily disconnected. In case of emergency, the disconnecting device is the mains lead. This will either be the mains lead connected to the product or, in the case of a cabinet-based system, the mains lead to the cabinet.
18. Only the correctly specified mains supply must be used.
19. Make sure the electrical cord is located so that it will not be subject to damage.
20. There are no user-serviceable parts in the camera. If the head is opened the warranty will be void. Only authorised service personnel may service this equipment.
21. Users must be authorised and trained personnel only; otherwise this may result in personal injury, and/or equipment damage and impaired system performance.



---

## Safety and Warning Symbols

The following are explanations of the symbols found on this product:



This product has been tested to the requirements of CAN/CSA-C22.2 No. 61010-1, 2nd edition, including Amendment 1, or a later version of the same standard incorporating the same level of testing requirements



The iXon Ultra and Life camera series require a Direct Current (DC) supply.



Refer to this guide before use.

## Manual Handling

Due to the delicate nature of some of the components within, care must be exercised when handling this product. Proper manual handling techniques are important when unpacking and installing the system to ensure that the integrity of the product is safeguarded and individuals involved are not exposed to unnecessary manual handling risks, such as:

- Lifting a load that is too heavy
- Poor posture or technique during lifting
- Dropping a load
- Lifting objects with sharp edges

## Shipping and Storage Conditions

Unpacking and Inspection:

- Carefully unpack the unit and retain packaging to return equipment for servicing.
- If the equipment appears damaged in any way, return it to sales outlet in its original packaging. No responsibility for damage arising from the use of non-approved packaging will be accepted.
- Ensure all items and accessories specified in **Section 1.5** are present.

If any items are missing, please contact your local sales representative.

---

## Section 1: Introduction

Thank you for choosing this Andor iXon camera. You are now in possession of a revolutionary new Electron Multiplying Charge Coupled Device (EMCCD), designed for the most challenging low-light imaging applications. This Hardware Guide contains useful information and advice to ensure you get the optimum performance from your new system.

---

## 1.1 Technical Support

If you have any questions regarding the use of this equipment, please contact the representative\* from whom your system was purchased, or:

### Europe

Andor Technology  
7 Millennium Way  
Springvale Business Park  
Belfast  
BT12 7AL  
Northern Ireland  
Tel. +44 (0) 28 9023 7126  
Fax. +44 (0) 28 9031 0792

### USA

Andor Technology  
300 Baker Avenue  
Suite # 150  
Concord  
MA 01742  
USA  
Tel. +1 (860) 290-9211  
Fax. +1 (860) 290-9566

### Asia-Pacific

Andor Technology (Japan)  
5F IS Building  
3-32-42 Higashi-Shinagawa  
Tokyo 140-0002  
Japan  
Tel. +81-(0)3-6732-8968  
Fax. +81-(0)3-6732-8939

### China

Andor Technology (China)  
Haitong Times Business Center,  
Building B2 West,  
No.11 West Third Ring North Road,  
Haidian District,  
Beijing,  
100089  
China  
Tel: +86 (0)10 5884 7900  
Fax. +86 (0)10 5884 7901

\* The latest contact details for your local representative can be found on the contact support page of our website [andor.oxinst.com/support/](http://andor.oxinst.com/support/)

---

## 1.2 Disclaimer

THE INFORMATION CONTAINED HEREIN IS PROVIDED "AS IS" WITHOUT WARRANTY, CONDITION OR REPRESENTATION OF ANY KIND, EITHER EXPRESS, IMPLIED, STATUTORY OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY WARRANTY OF MERCHANTABILITY, NON-INFRINGEMENT OR FITNESS FOR A PARTICULAR PURPOSE.

IN NO EVENT SHALL ANDOR BE LIABLE FOR ANY LOSS OR DAMAGE, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL, CONSEQUENTIAL OR OTHERWISE HOWSOEVER CAUSED WHETHER ARISING IN CONTRACT TORT OR OTHERWISE, ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE INFORMATION PROVIDED HEREIN.

PLEASE NOTE, AI TRANSLATIONS OF THIS USER MANUAL MAY RESULT IN INCORRECT INFORMATION. IN THE EVENT OF A CONFLICT WITH A FOREIGN LANGUAGE TRANSLATION, ANDOR'S ENGLISH TRANSLATION WILL PREVAIL.

## 1.3 Copyright and Protective Notices

The copyright in this document and the associated drawings are the property of Andor Technology Ltd. and all rights are reserved. This document and the associated drawings are issued on condition that they are not copied, reprinted or reproduced, nor their contents disclosed.

The publication of information in this documentation does not imply freedom from any patent or proprietary right of Andor Technology Ltd. or any third party.

## 1.4 Trademarks and Patent Information


iXon Ultra, Andor and the Andor logo are trademarks of Andor Technology Ltd. Andor Technology Ltd. is an Oxford Instruments company. All other marks are property of their owners.










### Manufacturers Information

**Andor Technology Ltd., Belfast, BT12 7AL, UK.**

## 1.5 Supplied Components

The Andor iXon Ultra and Life 897 are supplied with the following standard components:

Description		Quantity
	iXon Ultra 897 or Life 897 EMCCD Camera (The iXon Ultra 897 has an integrated C-mount shutter)	1

Description		Quantity	Description		Quantity
	USB cable	1 x 3 m		Mounting Posts (Ø1/2" x 80 mm long x 1/4-20 UNC)	2
	Multi I/O Timing cable ACZ-03463	1 x 3 m		C-mount cap (including O-ring)	1
	Power Supply Unit (PS-90)	1		iXon Ultra & Life 897 Quick Start Guide	1
	Country specific Power Cord	1		Performance Sheet	1
	Coolant hose connections (Coolant pipes)	2			

### 1.5.1 Optional Components

- Software SDK and/or Solis if ordered
- Andor Programmer guide to Andor Basic (if ordered)
- Optomask microscopy accessory, used to mask unwanted sensor area during Crop Mode acquisition (OPTOMASK).
- Re-circulator for enhanced cooling performance (XW-RECR)
- Oasis 160 Ultra compact chiller unit (ACC-XW-CHIL-160)
- C-mount to Nikon F-mount adapter (OA-CNAF)
- C-mount to Olympus adapter (OA-COFM)
- C-mount to T-mount adapter (OA-CTOT)

---

## 1.6 Camera Power Supply Unit (PSU)

The iXon Ultra and life require a Direct Current (DC) supply and are powered using an external 12 V PS-90 PSU. Please see "External Power Supply Requirements" on page 81 and "Camera Power Specifications" on page 81 for more information.

### 1.6.1 Working with Electronics

The computer equipment that is to be used with the camera should be fitted with appropriate surge/EMI/RFI protection on all power lines. Dedicated power lines or line isolation may be required for some extremely noisy sites. Appropriate static control procedures should be used during the installation of the system. Attention should be given to grounding. All cables should be fastened securely into place in order to provide a reliable connection and to prevent accidental disconnection.

The circuits used in the camera head are extremely sensitive to static electricity and radiated electromagnetic fields and should not be used (or stored close to) EMI/RFI generators, electrostatic field generators, electromagnetic or radioactive devices, or other similar sources of high energy fields. Types of equipment that can cause problems include Arc welders, Plasma sources, Pulsed-discharge optical sources, Radio frequency generators and X-ray instruments.

---

## 1.7 Prevention of Condensation

Condensation may form on the outside of the camera body if the temperature of the cooling water is too low or if the water flow rate is too high. The first signs of condensation will usually be visible around the connectors where the water tubes are attached. In such circumstances switch off the system, disconnect the power supply and carefully wipe the camera with a soft, dry cloth. It is likely there will already be condensation on the cooling block and cooling fins inside the camera. Please also carry out the following actions:

- Set the camera aside to dry for several hours before you attempt re-use
- Before re-use blow dry gas through the cooling slits on the side of the camera to remove any residual moisture
- Use warmer water or reduce the flow rate of water when you start using the device again
- Check Dew Point (refer to **Appendix B**)

## 1.8 EM Gain Ageing

It has been observed that some EMCCD sensors, more notably in cameras that incorporate L3Vision sensors from e2v, are susceptible to EM Gain fall-off over a period of time. This ageing effect applies to any EMCCD camera manufacturer that incorporates L3Vision sensors into their cameras - including the Andor iXon Ultra and Life 897 models.

A technical note entitled: 'EMCCD - RealGain™ & EMCAL™', which further explains this phenomenon, can be viewed on the Andor website: [andor.oxinst.com/learning/view/article/realgain.-anti-ageing-emcal](https://andor.oxinst.com/learning/view/article/realgain.-anti-ageing-emcal)

### 1.8.1 Minimizing EM Gain Ageing in your iXon Ultra/Life Camera

If left unchecked, EM Gain Ageing has the potential to significantly compromise the long-term quantitative reliability of EMCCD cameras. Andor has implemented innovative measures to stabilize the EM Gain on these sensors and ensure the long-term quantitative stability to the user. If these guidelines are followed EM Gain Ageing can be minimized and should not present any real problem to the user.

More details of this ageing effect and Andor's solutions can be found on **Section 4.1.7**. Some of the guidelines to minimize the EM Gain ageing process are listed below:

- **Do not use EM Gain values greater than is necessary to overcome the read noise.** A gain of x4 or x5 the rms read noise (accessible from the spec sheet or performance sheet) is more than sufficient to render this noise source negligible. In practice, this can be achieved with EM Gain of less than x300 at 10 MHz and x450 for 17 MHz operation. Pushing gain beyond this value would give little or no extra Signal to Noise benefit and would only reduce dynamic range.
- **Only select the extended EM Gain scale of x1000 for single photon counting applications** and always ensure that the signal falling onto the sensor is within the regime of low numbers of photons per pixel.
- Turn down the gain when the camera is not acquiring.
- Try not to over-saturate the EMCCD sensor.

---

## 1.9 Minimizing Particulate Contamination

It is important that particulate contamination of the exterior of the camera window is kept to a minimum so that images are kept free of 'shadowing' particles directly in the optical path. The iXon Ultra range comes equipped with an internal C-mount shutter- this is not present in the iXon Life range. Whilst not being required for frame transfer operation (which is a shutter-free readout mode) it is good practice to close the shutter when the camera is not in acquisition use for a reasonable period. It is also advisable to use the software to close the shutter when exposing the camera to the 'open environment' (i.e. removed from a microscope C-mount or focusing lens) whilst power is still flowing to the camera.

When exiting SOLIS the shutter (if fitted) will close automatically. We recommend that the C-mount opening of both the iXon Ultra and Life series is covered when the camera is not in use.

If there is evidence of particulate contamination on the front window it is possible to clean the window by blowing oil free dry air gently over the window surface. To ensure the shutter stays open (Ultra series only), unplug power from the camera when Solis is running, and the shutter is open. Note that exiting Solis abnormally may leave the shutter in the open state.

## 1.10 Software

The iXon Ultra and Life series can be supplied with Andor [Solis](#), [iQ](#) or [SDK](#) software. They are also compatible with a range of [third party software](#) options that support optimized acquisition control and analysis functionality.



---

## Section 2: Product Overview

### 2.1 iXon Ultra and Life 897

The iXon Ultra and Life 897 combine a 512 x 512 EMCCD sensor, USB 2.0 and optimized electronics of the iXon platform to provide the highest performance and reliability available for an EMCCD camera.




Figure 1: iXon Ultra 897 EMCCD Camera

Features of the iXon Ultra and Life 897 models include:

- 512 x 512 active pixels / 16  $\mu\text{m}$  pixel size
- 17 MHz readout - 26 fps
- 56 fps @ 512 x 512 (full frame), 569 fps @ 128 x 128 (ROI, Optically Centred Crop Mode)
- Single photon sensitivity
- Backilluminated > 90%  $\text{QE}_{\text{max}}$
- USB 2.0 Connectivity
- Ultravac™ Technology
- Cooling to -100°C iXon Ultra 897, (iXon Life 897 -80°C)

## 2.2 Power and Signal Connections

The power and signal connections are located on the base plate of the iXon Ultra and Life 897 models:




Figure 2: Power and Signal Connections of the iXon Ultra and Life 897 models

**Power connector:** For connection to the Power Supply Unit (PSU). Note the connector is keyed and has a locking action.

**USB 2.0:** A USB connection should be made to the PC with the supplied USB 2.0 cable – this type (with cable ferrite) was used during EMC testing. An alternative cable with locking adaptor is available.

**Ext I/O:** The iXon Ultra and Life models are supplied with an ACZ-03463 cable for the external I/O connector. This provides industry-standard output BNC connectors:

- Fire (please refer to Section 5.2)
- Shutter (see Section 5.3)
- Arm (please refer to Section 5.2)
- Ext. Trig (External Trigger Input) (please refer to Section 5.2)

These are used to send and receive Trigger and Fire signals. The outputs (Fire & Shutter) are CMOS compatible & series terminated at source (i.e. in the camera head) for a 50  $\Omega$  cable.

### NOTES:

1. The cable termination at the customer end should be high impedance ( $>1\text{ K}\Omega$ ) as an incorrect impedance match could cause errors with timing and triggering.
2. The External Trigger Input is TTL level, CMOS compatible and has 470  $\Omega$  impedance.
3. Signal diagrams of these connections are shown on Section 2.5. The interfaces and internal circuits of the iXon Ultra and Life models are rated as SELV (Safety Extra Low Voltage). All interfacing equipment should use SELV voltage and current levels.
4. OutputDAC1 and OutputDAC2 (iXon Ultra only) are 16 bit DAC outputs that can be configured by the user to be up to approximately 10.1 Volts. Maximum output current that can be drawn is 10 mA. Note DAC functionality is no longer available on models supplied after April 2022.
5. +5 V Output is a 5 V supply to signal to the user that the camera is powered up. The maximum current that can be drawn from this is 500 mA

6. I/O bits (8 off) (iXon Ultra only) are user programmable and can either be inputs or outputs. When being used as inputs these default to being weakly pulled high. The maximum low level input voltage is 1.5 V and the minimum high level input voltage is 3.5 V. As outputs the maximum "high" level output current that can be drawn is 0.03 mA and the maximum "low" level current that each output can sink is 10 mA.
7. I<sup>2</sup>C: I<sup>2</sup>C connection point- 2 x16 bit DACs, and 8 digital i/os available on the External I/O (iXon Ultra only). To access these connections requires an advanced cable (ACZ-03464) to connect to the 26 way High density D connector.

## 2.3 Camera Link (iXon Ultra only)

The iXon Ultra is equipped with a Base Configuration (3-tap interface) Camera Link output which conforms to the specification defined by the Automated Imaging Association (AIA). This provides access to the camera data output with very low latency. Note that this is an OUTPUT ONLY e.g for use with a Camera Link frame grabber or custom embedded applications.

A Technical Note "An Overview of the Camera Link Output Function for the iXon Ultra" is available at: [andor.oxinst.com/learning/view/article/camera-link-output](https://oxinst.com/learning/view/article/camera-link-output)

## 2.4 External I/O

The iXon Ultra and Life models have a 26 way high density D-type connector:




Figure 3: I/O connector (26-Way D-type 164A21019X)

Table 1: Pinouts for D-Type connector (on ACZ-03463 cable).

1	External Trigger	10	I/O data bit 0*	19	5 V Out
2	Trigger Invert	11	I/O data bit 1*	20	GND
3	GND	12	I/O data bit 2*	21	I <sup>2</sup> C Data
4	Output DAC 1*	13	I/O data bit 3*	22	I <sup>2</sup> C Clock
5	Output DAC 2*	14	I/O data bit 4*	23	Shutter Control Output
6	GND	15	I/O data bit 5*	24	Arm Output
7	Frame Output	16	I/O data bit 6*	25	GND
8	Fire Output	17	I/O data bit 7*	26	GND
9	Reserved Output	18	GND		

\* iXon Ultra models only. Note DAC functionality is no longer available on models supplied after April 2022.  
Note: Pin 9 is a reserved pin as interline shift is not used on the iXon.

## 2.4.1 Additional Cables (iXon Ultra Only)

The iXon Ultra is provided with a cable ACZ-03463 for the external I/O connector as standard. This provides industry-standard output BNC connectors wired to the Fire, External Trigger, Shutter and Arm. There are also two optional cables available for the iXon Ultra, which provide additional connectivity (ACZ-03453 and ACZ-03454).

The I/O timing cable (ACZ-03453) has the same BNC connections (Fire, External Trigger, Shutter and Arm) that the standard cable provides. However, the cable also has an additional 5 pin Fischer connector which links to the I<sup>2</sup>C data and I<sup>2</sup>C clock pins.

Whereas the optional I/O timing cable (ACZ-03454) has a total of 17 BNC connections. It includes the Fire, External Trigger, Shutter and Arm connections that the standard cable provides as well as the 13 additional connections highlighted in Table 2 below. However, it should be noted that no connections to the I<sup>2</sup>C data and I<sup>2</sup>C clock pins are included.

Table 2: Pinouts for D-Type connector (on ACZ-03454 cable).

1	External Trigger	10	I/O data bit 0	19	5 V Out
2	Trigger Invert	11	I/O data bit 1	20	GND
3	GND	12	I/O data bit 2	21	I <sup>2</sup> C Data
4	Output DAC 1*	13	I/O data bit 3	22	I <sup>2</sup> C Clock
5	Output DAC 2*	14	I/O data bit 4	23	Shutter Control Output
6	GND	15	I/O data bit 5	24	Arm Output
7	Frame Output	16	I/O data bit 6	25	GND
8	Fire Output	17	I/O data bit 7	26	GND
9	Reserved Output	18	GND		

\* iXon Ultra models only. Note DAC functionality is no longer available on models supplied after April 2022.

## 2.5 Signal Diagrams

### 2.5.1 iXon Ultra and Life Input & Output Timing Hardware




Figure 4: External Trigger

### 2.5.2 External Trigger Input (at connector)

- $V_{IH}$  (High level input voltage, minimum) = 2.2 V
- $V_{IL}$  (Low level input voltage, maximum) = 0.88 V
- Change to falling edge trigger by connecting "Trigger Invert Input" pin to Ground with  $\leq 500 \Omega$




Figure 5: Fire & Shutter

#### Fire, Shutter, Arm & Frame Output (at connector)

- Use 50 Ω cable + high impedance input ( $>1 \text{ k}\Omega$ )
- Drive both TTL & CMOS unterminated input

---

## Section 3: Installing the iXon Life and Ultra 897

### 3.1 PC Requirements

Install the camera software before first connecting the camera – this will ensure that USB drivers are available when required.

There are no restrictions on the order in which components are connected. It is best to allow a few seconds from camera power on (using either the button or a mains switch) to starting Solis in order for the camera to be recognised by the PC.

- 3 GHz Quad Core or 2.4 GHz multi core processor
- 2 GB RAM
- 100 MB free hard disc to install software (at least 1 GB recommended for data spooling)
- USB 2.0 High Speed Host Controller capable of sustained rate of 40 MB/s
- Windows 10, 11 or Linux

### 3.2 Connecting the Camera

1. Attach the camera to lens or optical system using the camera C-mount interface as required
2. Insert the 12 V DC power cable from the power supply into the power connector on the bottom plate of the camera, ensure the orientation is correct. NEVER forcibly insert the connector.
3. Connect the supplied USB 2.0 cable between the USB 2.0 connector on the camera and the corresponding slot on the PC. (A USB socket on the rear of a desktop machine is preferred to ensure full performance). Only use USB 2.0 cables supplied by Andor as performance cannot be assured with other models.
4. Switch the camera ON using the ON/OFF switch. You should hear an audible confirmation (camera start-up tone).

**Note:** The iXon Ultra and Life models have a power switch on the camera head for convenience.

5. The supplied Multi i/o cable may be required depending on the measurement being carried out. Refer to **Section 2.4** for details.
6. The camera can achieve stated performance with air cooling using the internal fan – Water cooling is also available see **Section 3.3** for details.

### 3.3 Cooling the CCD

Heat is generated by the sensor during normal operation which if not addressed may have a significant adverse effect on performance (e.g. signal to noise ratio and sensitivity) due to increased dark current noise. The iXon range make use of a multi-stage (3: iXon Life; 4 iXon Ultra) Peltier cooling assembly (thermoelectric cooler, TEC), which utilizes the thermoelectric effect to rapidly cool the sensor down to the stable operating temperature. A TEC has a cold side (in contact with the sensor) and a hot side. Temperature control components regulate the cooling of the camera and ensure that a stable temperature is maintained between and throughout measurements.

The iXon Ultra and Life models can use either forced air cooling- using the in-built fan, or water cooling for enhanced cooling performance (refer to **Section 3.4** for connection information). When using water cooling, a recirculator or a chiller can be purchased from Andor to provide a convenient and effective heat dissipation.

---

### 3.3.1 Sources of Heat Generation

In normal operation, clocking the image and storage regions of the EMCCD sensor, along with clocking the register, generates heat. The resistive heating process is dependent on the amplitude and frequency of the clocks, therefore the faster a sensor is clocked, the more heat is generated. A TEC is capable of providing a temperature difference (delta) between its cold and hot side that is dependent on the wattage of heat at the cold side (where the CCD is located). Therefore, the minimum temperature the CCD can achieve is dependent on the heat produced by the sensor.

Equally, if the hot side of the TEC can be maintained at a lower temperature, then the cold side will also be at a lower temperature, as the same delta is maintained. Air cooling the hot side achieves this by blowing air over the camera heat sink. The limitation of this is that the level of cooling is dependent on the temperature of the ambient environment. In the majority of cases this is sufficient, but if deeper cooling is necessary a water chiller or recirculator can be employed. This will keep the hot side of the TEC at a lower temperature as the heat can be transmitted to the water more efficiently and the water temperature can also be controlled.

## 3.4 Connecting a Cooling System

### 3.4.1 Hose Connections

Two barbed coolant hose inserts are supplied as standard with the iXon Ultra and Life camera series, suitable for connection to 6 mm (0.25") internal diameter soft PVC tubing / hose.

**Recommended tubing:** 10 mm (0.4") outside diameter, i.e. a wall thickness of 2 mm (0.08").

Alternative hose dimensions and materials should be thoroughly tested to ensure a leak tight seal is achieved with the barbed inserts.

### 3.4.2 Coolant Recommendations

Is recommended that de-ionized water (without additives) is used as the coolant to prevent deposits forming. Some mains supply water is heavily mineralized (i.e. "Hard") which could cause deposits in the water circuit inside the camera. This can reduce the flowrate and cooling efficiency.

The specified cooling performance of the camera can be achieved with coolant flow rates of >0.75 litres per minute, the maximum recommended pressure of coolant circulating through the camera head is 2 bar (30 PSI). In the event that replacement hose inserts / barbs are required, please contact your local Andor representative. **CAUTION:** Always ensure that the temperature of the liquid coolant circulated through the camera head is above the dew point of the camera ambient. Use of coolant at or below the dew point will result in permanent damage to the camera head, due to formation of condensation on internal components (refer to Appendix B).

---

### 3.4.3 Connecting the Coolant Hoses

1. Secure the hose to the barbs on the hose inserts
2. Click both hose inserts into the quick-release couplings on the side of the camera
3. Ensure the hose inserts are located securely in place



Figure 6: Hose inserts and quick release coupling

### 3.4.4 Removing the Coolant Hoses

**CAUTION:** Before attempting to remove the hose connections, ensure that all water has been drained from the hoses and the coolant channel within the camera head. Care must be taken to avoid permanent damage to the camera system resulting from either leakage of coolant during connection / removal of hoses or spillage of any residual coolant contained within the camera head once the hoses have been removed.

1. Depress the collar on the quick-release couplings
2. The hose insert connections will now be released and can be disconnected from the quick-release coupling



---

## 3.5 Mounting Posts

There are 4 pairs of mounting post positions on all four sides of the camera. These can be used to mount the camera if the C-Mount is not used, or to mount accessories. Each pair of holes has a 2.0" spacing.

Note: A bag containing two Ø1/2" x 80 mm long x 1/4-20 UNC posts is included with all kits.




Figure 7: Attaching the Mounting Posts

### 3.5.1 Attaching Mounting Posts to the Camera

1. Carefully remove the black grommet(s) as shown in Figure 7.
2. Screw each mounting post into the exposed mounting hole
3. Tighten using a screwdriver shank through the hole in the mounting post

**Note:** Store the blanking grommets so they may re-installed if the mounting posts are not in use.

---

## 3.6 Installing Andor Solis Software for Windows 10

1. Terminate & exit any applications which are running on the PC.
2. Insert the Andor Solis CD. The InstallShield Wizard should now start. If it does not start automatically, run the setup.exe file directly from the CD.
3. Select appropriate location for installation of software and drivers on your computer / network.
4. If prompted, select iXon Ultra/iXon Life.
5. Continue installation and restart your computer - when prompted - to successfully complete the installation.
6. The shortcut icon for Solis will appear on the desktop on re-start.
7. The iXon Camera is now ready to be connected to a PC / laptop and powered on.

## 3.7 New Hardware Wizard

When the iXon camera is connected to a PC for the first time, the **New Hardware Wizard** screen will appear.

1. Select the 'No, not this time only' option then click **Next>**.
2. Select the 'Install from a list or specified location (Advanced)' option then click **Next>**.
3. Navigate to the directory where the Andor Solis software was installed to on the PC, then click **Next>** so that the **Installation Wizard** can start.
4. Click the **Finish** button to complete the installation.

**Note:** If the camera is connected to a different USB port, steps 1 – 4 will have to be repeated on the first connection only.

5. A system message will appear to indicate that the device has been successfully installed.

**Note:** You can check that the iXon camera is correctly recognized and installed by opening the **Device Manager (Devices and printers)** in **Windows, Control Panel**. The iXon camera will show under the **Devices** list.

**Note:** On the first startup of Solis, you may be required to direct the software to the iXon Ultra/iXon Life drivers. If so, select the directory that Andor Solis was installed to.

---

## 3.8 Start-up Dialog

On start-up of Solis software a dialog may appear (similar to that shown below) if multiple cameras are connected to your PC.




Figure 8: Start Up Dialog Menu

1. Highlight the Andor iXon Ultra/ iXon Life camera (The Serial Number can be found on the label on the camera)
2. Click OK to continue with the selected camera.

## Section 4: Features and Functionality

### 4.1 EMCCD Operation

#### 4.1.1 Structure of an EMCCD

Advances in sensor technology have led to the development of a new generation of ultra-sensitive, low-light Electron Multiplying Charged Coupled Devices (EMCCDs). At the heart of your camera is the latest EMCCD, a revolutionary technology, capable of single photon detection. An EMCCD is a silicon-based semiconductor chip bearing a two-dimensional matrix of photo-sensors or pixels. This matrix is usually referred to as the image area. The pixels are often described as being arranged in rows and columns, the rows running horizontally and the columns vertically. The EMCCD in the camera is identical in structure to a conventional Charged Coupled Device (CCD) but with the shift register extended to include an additional section, the Multiplication or Gain Register as shown in Figure 9 below:




Figure 9: Structure of the EMCCD Sensor.

During an acquisition using a conventional Frame Transfer CCD (FT CCD), the image area is exposed to light and an image is captured. This image, in the form of an electronic charge, is then automatically shifted downwards behind the masked region of the chip before being read out. To read out the sensor, charge is moved vertically into the readout register, and then horizontally from the readout register into the output node of the amplifier. As stated previously, the readout register is extended to include the **multiplication (gain) register**. The amplification occurs in this register through the scheme highlighted in Figure 10 below. When moving charge through a register, there is a very tiny, but finite probability that the charges being transferred can create additional charge by a process known as “impact ionization”. Impact ionization occurs when a charge is clocked and has sufficient energy to create another electron-hole pair in the conduction band and by this process amplification occurs. To make this process viable, EMCCD’s optimize the process in two ways:

1. The probability of any one charge creating a secondary electron is increased by giving the initial electron charge more energy. This is typically done by replacing one of the electrodes (phases) of this readout section with two electrodes. The first is held at a fixed potential, and the second is operated as normal, except that much higher voltages are used than are necessary for charge transfer alone. The large electric field generated between the fixed voltage electrode and the clocked electrode is sufficiently high for the electrons to cause “impact ionization” as they transfer. This impact ionization causes the generation of new electrons, i.e. multiplication or gain.
2. The EMCCD is designed with hundreds of cells or pixels in which impact ionization can occur and although the probability of amplification or multiplication in any one pixel is small (only around  $\times 1.01$  to  $\times 1.015$  times) over the entire length of the EM register, the probability is very high and substantial gains of up to thousands can be achieved.




Figure 10: Gain register operation

---

## 4.1.2 EM Gain & Read Noise

As explained in 4.1.1, EMCCD sensors allow the detected signal to be amplified on the actual sensor itself before being readout through the output amplifier and digitized by the Analog to Digital (A/D) converter. The reason that this on-chip-multiplication process gives such a significant improvement in low-light detection is that it negates the effect of any electronic noise that may be generated by the readout electronics.

All CCD cameras have an associated minimum electronic noise floor, which is often termed the **Read Noise** of the system. Read noise is produced during the readout process mostly by the output amplifier, but also has contributions from the digitization electronics. This sets the minimum signal level that can be detected by the camera, as any signal level below the read noise level will be indistinguishable from the read noise itself.

Read noise has therefore been the major limiting factor for low-light level detection in CCDs for many years until the introduction of EMCCD cameras by Andor Technology in 2000. By applying EM Gain, a weak signal that would otherwise be indistinguishable from the read noise can be amplified above the read noise level and thus be read out as a useful signal. This amplification of the signal before being read out effectively reduces the read noise level of the camera, and even at relatively modest EM Gain settings the effective read noise can be reduced to less than 1 electron rms.

**NOTE:** Read noise increases with increased readout rate, therefore the application of EM Gain provides an advantage at higher readout rates, as any increase in the read noise can be overcome simply by increasing the EM Gain. For example, an iXon Ultra 897 typically has a read noise of 50 electrons rms when reading out at 10MHz. This can easily be reduced to < 1 electron by applying > x50 EM Gain.

### 4.1.3 EM Gain ON vs EM Gain OFF

Signal to Noise (S/N) plots derived from the specifications of the back-illuminated iXon Ultra EMCCDs (Figure 11), read out at 10MHz for a photon wavelength at which the Quantum Efficiency (QE) of the sensor is assumed to be 90%. Such plots are very useful to gauge at what signal intensity it becomes appropriate to use EM Gain to increase S/N.

It is clear that at 10MHz readout, one needs relatively intense signals of > 2900 photons / pixel before it becomes advantageous to operate with EM Gain off. Note that the "ideal" curve represents a pure Signal to Shot Noise ratio and is shown for reference – if the camera had no sources of noise, this is what the curve would appear like. Even with EM Gain turned on, we encounter uniformly lower signal to noise than the ideal curve. This is due to the influence of Multiplicative Noise, which has the effect of increasing the shot noise by a factor of  $\sqrt{2}$  or  $\sim 1.41$ .




Figure 11: Signal to noise plots EM Gain ON vs. EM Gain OFF for backilluminated iXon Ultra EMCCDs at 10 MHz readout speed

S/N plots derived from the specifications of the back-illuminated iXon Ultra EMCCDs at 1 MHz (slower frame rate operation) read out either with EM Gain ON or, alternatively, for the iXon Ultra series only, through the conventional amplifier (i.e. standard CCD operation) are shown in Figure 12. Again, this plot assumes a photon wavelength at which the QE of the sensor is 90%. This figure applies specifically to models where the user has the choice of either EMCCD or conventional amplifiers.

At these slower speed operations, when one has the choice to read out as a "conventional" CCD, it can often be advantageous to do so in order to achieve improved signal to noise. The plots below show that the cross-over point is at  $\sim 42$  photons/pixel, at which it is still advisable to read out through the EM amplifier with Gain applied.




Figure 12: EM Gain ON vs. Conventional Amplifier signal to noise plots for back-illuminated iXon Ultra EMCCDs at 1 MHz readout speed.

## 4.1.4 Multiplicative Noise Factor and Photon Counting

It is impossible to determine the exact gain a detected signal charge traversing the EM Gain register will acquire, due to the stochastic nature of the processes which produce EM Gain. However, it is possible to calculate the probability distribution function of output charge for a given input charge.

At high gain levels ( $> \times 30$ ) this uncertainty introduces an additional noise component called **Multiplicative Noise**. This noise source is only present in signal amplifying technologies and is a measure of the uncertainty inherent in the signal multiplying process. For example, during each transfer of electrons from element to element along the gain register of the EMCCD, only a small probability exists that the process of impact ionization will produce an extra electron during that step. This is a small probability, but when executed over  $> 590$  steps it results in a very large overall EM Gain. However, the downside to this process results from the probabilities. Due to this, there is a statistical variation in the overall number of electrons generated by the gain register from an initial charge packet. This uncertainty is quantified by a parameter called "**Noise Factor**" and detailed theoretical and measured analysis has placed this Noise Factor at a value of  $\sqrt{2}$  (or 1.41:1) for EMCCD technology.

**Note:** This noise source is significantly greater for the Multi Channel Plate (MCP) of ICCDs than for the gain register of EMCCD. ICCDs have noise factors typically ranging from 1.5 to  $>2$ .

Multiplicative noise is an additional form of noise that must be taken into account when calculating Signal/Noise for these types of detectors. However, one way to better understand the potential effects of this noise source is in terms of an addition to the shot noise of the system. Extra multiplicative noise has the same form as shot noise as each noise type results in an increase in the variation of number of electrons that are read out of the sensor (under constant uniform illumination).

Multiplicative noise can be thought to contribute directly to the overall shot noise, in that one should multiply the **Shot Noise** by the **Noise Factor** when calculating overall noise. Simply put, multiplicative noise does not in any way reduce the average signal intensity or reduce the number of photons that are detected. It simply increases the degree of variation of the signal around the mean value, in addition to the variation that already exists from the shot noise (variation from pixel to pixel or from frame to frame). This additional variation to the signal intensity is represented in Figure 13 as a signal intensity profile.




Figure 13: Signal Intensity Profile

Within the framework of less than 1 electron falling on a pixel in a single exposure, the EMCCD can be used in **Photon Counting Mode**. In this mode, a threshold is set above the ordinary amplifier readout and all events are counted as single photons. In this mode, with a suitably high gain, a high percentage ( $>90\%$ ) of the incident photons can be counted without being affected by the Noise Factor effect.



### 4.1.5 EM Gain Dependence and Stability

EM Gain is a function of the EM voltage and of the sensor operating temperature. When the user applies gain through the software, it is the EM voltage in the gain register that is varied. As can be seen from **Figure 14** below, the dependence of EM Gain on EM voltage is sharp (note the logarithmic scale). This arises because the signal electrons acquire energy as they are accelerated through the EM electric field, and once this field strength reaches the threshold needed to overcome the bandgap energy, the impact ionization rate rises rapidly. Historically, this sharp dependence has meant that the software control of EM Gain in all EMCCD cameras has been via a non-linear scale with most of the amplification occurring within a relatively small portion at the top of the overall scale. Thus, considerable fine tuning by the user to determine an optimal gain setting has been required, and even then the actual gain is determined only through measurement of a stable light source, with and without gain applied.




Figure 14: EM Gain vs EM clock voltage

The effect of sensor temperature on EM Gain is shown in **Figure 15**. This dependence arises primarily from photon scattering of electrons when they are accelerating in the EM electric field. The scattering causes a loss of energy, which increases with temperature. To make up this loss and maintain EM Gain, a larger EM electric field must be used at higher temperatures. As can be seen from **Figure 15**, EM Gains well in excess of x1000 can be achieved at low temperatures. However, it is not recommended to use gains above x1000 because such high gains can cause significant ageing of the gain register (see EM Gain Ageing in **Section 4.1.7**).




Figure 15: EM Gain vs sensor cooling temperature

## 4.1.6 RealGain™: Real and Linear gain

Andor has successfully converted the relationship between EM Gain and the EM clock voltage setting into a linear one through a detailed analysis of the complex EM voltage dependence. The actual EM Gain can therefore be selected directly from a linear scale displayed in software.

In effect one can select the best gain to overcome noise and maximize dynamic range. Also, although EM Gain is temperature dependent, Andor's linear and real gain calibration extends to any EMCCD cooling temperature. Selecting x300 EM Gain at -50°C, or at -100°C gives the same x300 actual gain.




Figure 16: RealGain™ calibration in the iXon Ultra – the same linear relationship holds across all cooling temperatures

---

## 4.1.7 EM Gain Ageing: What causes it and how is it countered?

As already noted in **Section 4.1.5**, EMCCD sensors can suffer from EM Gain ageing. This is the phenomenon whereby the EM Gain falls off over a period of time, when operating at the same clock voltage and cooling temperature. This ageing effect appears to be dependent on the amount of charge that is passed through the gain register, combined with the actual EM electric field strength that it is transferred through. It seems to be very strongly dependent on the EM electric field strength. Therefore, when operating at high EM Gains the ageing rate can be disproportionately greater. Fortunately, it has been observed that this ageing effect itself decreases with time, meaning that, with proper use, the device should remain useful for many years. As part of Andor's EMCCD production process, all sensors are conditioned so much of the "shorter-term ageing" has already occurred prior to calibration and setting of the EM Gain.

The rationale for this ageing effect is not fully understood, but it is assumed that accelerating charge through the high electric fields is causing a tiny fraction of that charge to become permanently embedded in the insulator (typically silicon dioxide) between the EM electrode and the active silicon. This slow build-up of charge effectively reduces the field strength produced by the electrode. The signal electrons, therefore, experience a lower accelerating potential, which subsequently produces fewer secondary electrons from the impact ionisation process resulting in less electron multiplication and, in effect, a lower EM Gain.

In order to minimise the effect of EM Gain ageing it is recommended that the following guidelines are always adhered to:

- Do not use EM Gains greater than necessary to overcome the read noise (please refer to **Figure 15** and **Figure 17**). A gain of x4 or x5 of the root-mean-square read noise (accessible from the performance sheet) is normally more than sufficient to render this noise source negligible. In practice, this can always be achieved with EM Gain of less than x500 (much less for the slower readout speeds). Pushing gain beyond this value would give little or no extra S/N benefit and would only reduce dynamic range.
- Only select the extended EM Gain scale of x1000 when **single photon counting** and always ensure that the signal falling onto the sensor is within the regime of low numbers of photons per pixel.
- Turn the EM Gain OFF when not in use.
- Try not to over-saturate the EMCCD detector.

Andor's linear and quantitative gain calibration scale (**RealGain™**) is one of a number of systems that prevent conditions from being selected that may result in accelerated sensor ageing in addition to allowing the optimum imaging performance to easily be achieved.

## 4.1.8 Gain and Signal Restrictions

Part of the measures taken to prevent premature ageing of the sensor has been to use temperature compensated real gain limits, coupled with signal intensity feedback (after EM amplification). This ensures that the user is unable to apply excessive gain and/or signal, any more than is necessary to render the read noise floor negligible for a given signal intensity and readout speed. Secondly, when not actually acquiring data, as for example during "keep clean" cycles or when outside a selected sub-image area, Andor EMCCDs have been internally configured to prevent any unwanted signal entering the EM Gain register. Together, these measures ensure that the rate of EM Gain ageing is significantly reduced.

## 4.1.9 EMCAL™

Andor has developed, in the iXon Ultra and Life models, a unique and patented method of user-initiated EM Gain self-recalibration - EMCAL™.

The reduction in EM Gain over time can be corrected by using the EMCAL™ self-recalibration process, which is very easily initiated by the user. Please check the Andor website for the latest EMCAL™ routine.

This process uses the in-built temperature compensated linear gain scales to reset the EM Gain calibration to reflect the true values requested on the software scale, giving the RealGain™ values. This means optimal signal to noise ratio, maximum dynamic range and prolonged system longevity (See Figure 17).




Figure 17: Ageing profile of an Andor backlit EMCCD. Test conditions: 24/7 operation; 30 frames/sec; x90,000 electrons per pixel through gain register; ~ 200,000 pixels illuminated

---

## 4.2 Sensor Cooling Considerations

### 4.2.1 Thermoelectric Cooling

The iXon Ultra and Life models use a Thermoelectric (TE) cooler, which alongside the Ultravac™ (see [Section 4.3](#)) design allows a minimum cooling temperature performance unrivalled by other manufacturers. TE coolers are small, electrically powered devices with no moving parts, making them reliable and convenient. A TE cooler is actually a heat pump, i.e. it achieves a temperature difference by transferring heat from its “cold side” (the EMCCD sensor) to its “hot side” (the built-in heat sink). Therefore, the minimum absolute operating temperature of the EMCCD depends on the temperature of the heat sink. The maximum temperature difference that the iXon Ultra and Life TE device can attain is dependent on the following two factors:

- Heat generated in the sensor will vary with readout mode
- TE-cooler “hot side” Temperature

For every 7 degrees cooling, dark current approximately halves. However, for a given temperature the actual dark current can vary by more than an order of magnitude from device to device. Cameras are specified in terms of the minimum dark current achievable, rather than minimum temperature.

### 4.2.2 Readout Mode and Heat Generation in the EMCCD

Heat is generated by the EMCCD sensor, primarily from the Vertical Clocks. These can be either the ‘storage’ or ‘image’ clocks and are used to transfer charges in the EMCCD vertically, down to the ‘shift register’ to be amplified and digitized. Increasing the horizontal readout rate from 10 MHz to 17 MHz has no impact on these, however the reduced time performing horizontal readout means that a higher proportion of a readout cycle is spent using vertical shifts i.e. more heat generated in the sensor per unit time.

Boosting the amplitude of the vertical clocks will generate the most heat, however any mode which increases the number of vertical clocks per unit time may also impact cooling. Readout settings which may affect minimum cooling are:

- Vertical shift speed
- Vertical binning e.g. region of interest or Full Vertical Binning (FVB)

Moving from Full image readout to FVB will reduce minimum achievable EMCCD temperature by approximately 10°C.

**Note:** The OptAcquire “Fastest Frame Rate” setting in Solis uses the fastest vertical shift speed and maximum vertical clock amplitude boost. This limits the minimum EMCCD cooling achievable.

### 4.2.3 Heatsink “hot side” Temperature

The minimum sensor cooling depends on the temperature of the heat sink. The maximum EMCCD cooling will be achieved when the “hot side” is the coldest possible. This is usually achieved with water cooling. Water cooling maintains the “hot side” at the water temperature irrespective of ambient air temperature.

## 4.2.4 Cooling Options

The heat that builds up on the heat sink must be removed and this can be achieved in one of the two following ways:

1. **Air cooling:** a small built-in fan forces air over the heat sink
2. **Water cooling:** external water is circulated through the heat sink using the water connectors on the head and this can take one of the following forms:
  - Recirculation
  - Chilling

All Andor iXon Ultra and Life systems support both cooling options. Irrespective of which method is being employed, it is not desirable for the operating temperature of the EMCCD simply to be dependent on, or vary with, the heat sink temperature. Therefore, a temperature sensor on the EMCCD (combined with a feedback circuit that controls the operating current of the cooler) allows stabilization of the EMCCD to any desired temperature within the cooler operating range.

## 4.2.5 Adjusting the Sensor Operating Temperature

The Sensor Temperature can be adjusted in the Temperature Control Dialog, shown below, this is accessed via the **Hardware/Temperature** menu:




Figure 18: The Temperature Control Menu

Or by double clicking on the temperature display at the bottom left corner of the main SOLIS application window:

OFF

Once the cooler is switched ON the current temperature is displayed in the bottom left of the main application window and the background of the icon will turn red.

25°C

As the system cools, the temperature displayed will decrease in 1°C increments. Indication of sensor temperature stabilization in Solis software is given by a coloured indicator at the bottom left of the main window, changing from RED (target cooling temperature not stabilized) to BLUE (stabilized).

---

## 4.2.6 Cooling Fan Settings

The speed and operation of the cooling fan can also be software-controlled, which is a prerequisite when working with experimental configurations that are extremely sensitive to vibration. The vast majority of applications, including optical microscopy, can be used with the default highest fan speed, since the vibrations from the fan are minimal. However, some applications can be extremely sensitive to even the smallest of vibrations (such as when combining an optical set-up with patch clamp electrophysiology or atomic force microscopy) and it can be useful to either select a slower fan speed, or to temporarily turn off the fan altogether, for the duration of the acquisition. The cooling fan settings can be accessed through the **Hardware/FanControl** menu

As the fan is used to dissipate heat generated from the TE cooler to achieve EMCCD cooling, switching it off will require a higher EMCCD temperature to be chosen.

**Note:** It is recommended to achieve a stable cooling temperature after the fan has been switched off, otherwise transient dark current behaviour will be seen.

## 4.3 Sensor Readout Optimisation

To allow the camera to be optimized for the widest range of applications, it is important to have flexibility in the readout options available. Some of these include:

- Cooling (see Section 4.2)
- Sensor pre-amp settings (see Section 4.3.1)
- Variable horizontal readout rate (see Section 4.3.2)
- Variable vertical shift speed (see Section 4.3.3)
- Output amplifier selection (see Section 4.3.4)
- Baseline settings (see Section 4.3.5)
- Binning and Sub Image options (see Section 4.3.6)

### 4.3.1 Sensor Pre-amp Options

An EMCCD sensor can have a much larger dynamic range than can be faithfully reproduced with the Analogue/Digital converters and signal processing circuitry currently available. To overcome this shortcoming, and to access the range of signals from the smallest to the largest, as well as to optimize the camera performance, it is necessary to allow different pre-amplifier gain settings.

However, with regards to selecting something other than the default highest pre-amp (most sensitive) setting for applications, it is recommended that this only ever be carried out when faced with extremely challenging dynamic range concerns. It is very important, however, that for such high-dynamic range applications, the user applies even more care to the amount of EM Gain applied (high EM Gain can drastically reduce the true dynamic range of the camera). Ideally, for maximum dynamic range whilst maintaining improved Signal to Noise (S/N), the EM Gain setting should be set equal to the read noise at the readout speed selected (value obtainable from the performance sheet that comes with the delivered system).

Pre-amplifier gain selection in CCDs is traditionally used to trade off S/N vs dynamic range. A higher pre-amp setting means fewer electrons/count, resulting in a lower system noise floor, therefore a higher S/N. However, high pre-amp settings may not match well to the pixel well depth of the sensor, therefore a lower setting can be selected to meet the full well depth potential, e.g. a pre-amp setting yielding 1.5 e-/count may be selected to ensure that the 65536 digitization levels of a 16-bit A/D are closely matched to a 100,000 e- pixel well depth. A pre-amp setting of 1 e-/count, while giving a lower noise floor, would not harness the full 100,000 e- well depth within the 16 bit A/D.

---

The situation is not nearly as straightforward for EMCCDs because:

1. EM Gain overcomes readout noise and amplifies signals relative to the digitization noise (which is fixed for a given pre-amp setting).
2. Gain register pixels have a greater well depth than the imaging pixel well depth.

The latter point can be particularly confusing and has led to confusion in the field. What this has meant, is that we have set some of the lower pre-amp settings associated with the EM-output to match the extended well capacity of the gain register pixels (as reported by the sensor manufacturer e2v). This means that these pre-amp settings are designed to be used with EM Gain. Otherwise, the lower well capacity of the imaging pixels will saturate long before the A/D. This is why some users have been confused at not being able to reach the full ~65k counts of the 16-bit A/D channel, when they hadn't applied EM Gain.

Andor recommend using the default highest value pre-amp setting (e.g. Gain3 setting of the iXon Ultra 897E giving ~ 4 e-/count @ 10MHz) for most low-light applications. Most genuinely low-light applications are not limited by well capacity, as long as sensible EM Gain settings are applied (we recommend not exceeding x500 EM Gain, except for single photon counting experiments). Even at this highest pre-amp setting, the typical imaging pixel well depth will still be exceeded before the 16-bit A/D would saturate ( $180,000 \text{ e-} \div 4 \text{ e-/count} = 45,000 \text{ e-}$ ). The remainder of the 16-bit A/D range is therefore still free to be utilized by the extended well capacity of the gain register. For example, with an EM Gain of x300 (RealGain™), it would take 600 electrons in a pixel of the sensor to reach this A/D saturation limit. Say the QE is 80% at the wavelength of interest then this corresponds to maximum of 750 photons falling onto that pixel. That is perfectly satisfactory dynamic range for the vast majority of low-light imaging applications.

**NOTE: A Side effect of the new high speed ADC method chosen for the iXon Ultra, is that the full range of ADC codes 0- 65535 is not available. A margin has been added and the iXon will reach an ADC saturation that is lower than 65535. See Camera performance sheet for details.**

The core reason for recommendation of this pre-amp setting, even over the middle (Gain2) pre-amp setting, is that it implements an additional restriction as to how much charge is allowed to build up in the sensor. This in turn will help minimize the rate of EM Gain ageing (please see [Section 4.1.7](#) for further details on measures against gain ageing). However, some applications can be very demanding of dynamic range, and for those a lower pre-amp setting such as Gain2 is recommended. This ensures the A/D capacity is more closely matched to the well capacity of the gain register pixels and provides the maximum dynamic range. Also, as mentioned above, to maximize the true dynamic range of the camera tuning the RealGain™ gain setting to a value close to the value of the readout noise at the selected readout speed is recommended (e.g. if readout noise is ~ 50 electrons @ 10 MHz, set the EM Gain to x50 for maximum dynamic range).

## 4.3.2 Variable Horizontal Readout Rate

The **Horizontal Readout Rate** defines the rate at which pixels are read from the shift register. The faster the horizontal readout rate the higher the frame rate that can be achieved. The ability to change the pixel readout speed is important in order to achieve the maximum flexibility of camera operation, particularly in terms of dynamic range:

- **Slower readout rates:** Lower read noise and higher available dynamic range, but at the expense of slower frame rates.
- **Higher readout rates:** Highest frame rates, increase in read noise and reduced available dynamic range.

There are a number of different horizontal readout rates available on the iXon Ultra and Life models. Please refer to the performance sheet for readout rates available on your particular model.



---

### 4.3.3 Variable Vertical Shift Speed

The **vertical shift speed** is the time taken to vertically shift all pixels one row down, with the bottom row entering the shift register. The ability to vary the vertical shift speed is important for several reasons. It is possible, by using the different vertical speeds, to better synchronize the frame rates to external events such as a confocal spinning disc. Faster vertical shift speeds also have benefits such as lower Clock Induced Charge (CIC). A drawback with faster vertical shift speeds, is that the charge transfer efficiency is reduced, effectively reducing the pixel well depth. This is particularly important for bright signals, as a pixel with a large signal is likely to have some charge left behind if the vertical shift speed is too fast. This will result in degraded spatial resolution.

Slower vertical clocks ensure better charge transfer efficiency, thus giving maximum pixel well depth, but result in a slower maximum frame rate. To improve the transfer efficiency the clocking voltage can be increased using the vertical clock voltage amplitude setting. However, the higher the voltage, the higher the clock-induced charge. Thus, the user must make a measured judgement as to which setting works best for their situation.

- **For low CIC:** Use the fastest vertical shift speed that still transfers charge correctly (no image distortion), without having to select excess vertical shift voltage amplitude.
- **For maximum pixel well depth:** Use the slowest vertical shift speed, which will give an increase in CIC.
- **For maximum frame rate:** Use the fastest vertical shift speed and increase the vertical shift voltage amplitude to the minimum value that regains the full pixel well depth.
- **To reduce vertical smearing during very short exposure:** Use a faster vertical shift speed. This vertical smearing is due to the fact that light is still falling on the image area during the short time taken to transfer the charge from the image area into the storage area. If the actual exposure time is of a similar magnitude to this transfer time, then as pixels are shifted vertically through brighter regions of the image they will collect "extra" charge which will manifest itself as vertical streaking.

**Note:** For extremely short exposure times, a fast external shutter or pulsed light source may be required.

- **For short exposures (e.g. 1ms):** With high signal count and DC illumination, it may be necessary to increase the vertical clock voltage to ensure that the Keep Clean Cycle can fully remove the extremely high (saturated) signal that may have accumulated during the sensor readout phase.

## 4.3.4 Output Amplifier Selection

The iXon Ultra camera incorporates dual output amplifiers, an electron multiplying output amplifier and a conventional output amplifier. This increases the versatility of the camera:

- **EM amplifier:** Select for fast imaging in low-light conditions.
- **Conventional amplifier:** (iXon Ultra only) select where more light is available and a slower readout, with its associated lower read noise and higher dynamic range, is preferred.

The readout structure on sensors with both output amplifiers present is shown in Figure 19. From this it can be seen that, when reading out through the EM amplifier, accumulated charge will move to the right along the serial register and then into the EM Gain register. When the conventional output amplifier is selected, the charge to be read out will move along the serial register to the left then be transferred directly into the conventional output amplifier. This change in direction has the effect of producing mirror images when comparing raw data from the two output amplifiers. Some software packages will automatically reverse the image orientation of one of the output amplifiers to allow direct comparison of images. Consult your software manual to verify if this is the case.




Figure 19: Sensor readout structure

## 4.3.5 Baseline Optimization

### 4.3.5.1 Baseline Clamp

When acquiring data, small changes in the ambient temperature and/or in the heat generation of the driving electronics within the camera may cause some drift in the baseline level. This is most often observed during long kinetic series. The iXon Ultra and Life camera series employ a **Baseline Clamp** technique that holds the baseline to a predetermined level. Baseline Clamp corrects each individual image for any baseline drift by subtracting an average bias signal from each image pixel and then adding a fixed count level to ensure that the displayed signal level is always a positive number of counts.

**Note:** Baseline clamp is permanently on during normal operation.

### 4.3.6 Binning and Sub Image options

**Binning** is a process that allows charge from two or more pixels to be combined on the EMCCD-chip prior to readout. Summing charge on the EMCCD, and doing a single readout, gives better noise performance than reading out several pixels and then summing them in computer memory. This is because each movement of the charge through the readout amplifier contributes to the noise. There are two types of binning:

- **Vertical Binning:** Where charge from two or more rows of the EMCCD-chip are moved down into the shift register before the charge is read out. The number of rows shifted depends on the binning pattern selected. Thus, for each column of the EMCCD-chip, charge from two or more vertical elements is **summed** into the corresponding element of the shift register. The charge from each of the pixels in the shift register is then shifted horizontally to the output amplifier and read out.
- **Horizontal Binning:** Where charge from two or more pixels in the serial register are transferred into the output amplifier and read out as one combined data value. Thus, the charge from two or more of the horizontal elements is effectively summed into the output amplifier before being readout.

Combining both the vertical and horizontal binning methods produces "Superpixels". These consist of two or more individual pixels that are binned and read out as one large pixel. Thus, the whole EMCCD, or a selected sub-area, becomes a matrix of Superpixels, e.g.:




Figure 20: Representation of Superpixels

The horizontal and vertical binning parameters determine the dimensions of any superpixels created. On the one hand superpixels result in a loss of spatial resolution when compared to single pixel readout, but on the other hand they offer the advantage of summing data on-chip prior to readout thereby producing a better signal to noise ratio and a higher frame rate. All iXon Ultra and Life models offer completely flexible binning patterns which can be selected in the control software.

For the purpose of initial focusing and alignment of the camera, or to increase the readout speed, the user may wish to only readout a particular sub-area of the EMCCD to produce a Sub Image.

When a Sub Image has been defined, only data from the selected pixels will be digitized. Data from the remaining pixels will be discarded. The flexible configuration of the iXon Ultra allows the user to set the Sub Image area to any size and location on the EMCCD chip.

**Note:** Due to the wave shape presented by the sensor at fast clocking speeds, horizontal binning at 17 MHz does not give a linear relationship between 1x1 and other horizontal binning options (2x1, 4x1, etc.) Vertical binning is unaffected.




Figure 21: Vertical and Horizontal binning of two rows

- Step 1** Charge is built up on the sensor.
- Step 2** Charge in the frame is shifted vertically by one row, so that the bottom row of charge moves down into the shift register.
- Step 3** Charge in the frame is shifted vertically by a further row, so that the next row of charge moves down into the shift register, which now contains charge from two rows - i.e. the charge is vertically binned.
- Step 4** Charge in the shift register is moved horizontally (through the EM Gain register, if using the EM output amplifier) until the charge from the first data pixel is just about to enter the output node of the amplifier.
- Step 5** Charge in the shift / EM Gain register is moved horizontally by one pixel, so that charge on the end-most pixel of the shift register is transferred into the output node of the amplifier.
- Step 6** Charge in the shift register is again moved horizontally, so that the output node of the amplifier now contains charge from two pixels of the shift register - i.e. the charge has been horizontally binned.
- Step 7** The charge in the output node of the amplifier is passed to the analog-to-digital converter for each row and is read out.
- Step 8** Steps 5 - 7 are repeated until the shift register is empty. The process is repeated from Step 2 until the whole frame is read out.

## 4.4 Acquisition Options

### 4.4.1 Capture Sequence in Frame Transfer (FT) Mode

A number of acquisition modes are available to best suit your experimental demands. In **Frame Transfer (FT)** acquisition mode, the fastest performance can be delivered whilst maintaining optimal Signal to Noise. It achieves this through simultaneously acquiring an image onto the image area whilst reading out the previous image from the masked frame storage area. Thus, there is no time wasted during the readout and the camera operates with what is known as a 100% 'duty cycle'.




Figure 22: Capture sequence (FT mode)

- Step 1** Both Image and Storage areas of the EMCCD are fully cleaned out. This is known as a "Keep Clean Cycle" (KCC). Keep Clean Cycles occur continuously to ensure that the camera is always ready to start an acquisition when required. Further details of the Keep Clean Cycle are given later.
- Step 2** On receipt of a Start acquisition command the EMCCD stops the Keep Clean Cycle. This allows the image (photoelectric charge) to build up in the Image area of the EMCCD. The EMCCD remains in this state until the exposure time has elapsed, at which point the readout process starts.
- Step 3** The first phase of the readout process is to quickly shift the charge, built up in the Image area, into the Storage area. The time required to move the charge into the Storage area is approximately calculated as follows: (No. of rows in the Image area) x (vertical shift rate).
- Step 4** Once the Image area has been shifted into the Storage area the Image area stops vertically shifting and begins to accumulate charge again, i.e. the next exposure starts. While the Image area is accumulating charge the Storage area is being read out. This readout phase can take tens of milliseconds to seconds depending on the image size, readout pattern and readout speed.
- Step 5** On completion of the readout, the system will wait until the exposure time has elapsed before starting the next readout (i.e. returning to Step 3).

---

#### 4.4.1.1 Points to consider when using FT Mode

- In this mode, there are no Keep Clean Cycles between images during an accumulation or kinetics series as they are not necessary.
- This mode gives the fastest way to continually take images; however, the minimum exposure time is restricted to the time taken to read out the image from the Storage area.
- The accumulation cycle time and the kinetic cycle time are fully dependent on the exposure time and hence cannot be set via software.
- In external trigger mode there are no Keep Clean Cycles and the External trigger starts the "read out" phase. The exposure time is the time between external triggers and hence the user cannot set the exposure or cycle times. However, the user can define the amount of time between the external trigger event occurring and the readout starting. This can be useful in those situations where the TTL trigger occurs before the light event you are trying to capture. This effectively moves the exposure window in time, but the exposure time is still the period between trigger events.
- There is no need for a mechanical shutter. The exposure time (which is at least equal to the time to readout an entire image) is long compared to the time required to shift the image into the Storage area, and therefore image streaking will be insignificant.

## 4.4.2 Capture Sequence in Non-Frame Transfer Mode (NFT) with an FT EMCCD

It is also possible to operate an FT EMCCD in a Non-Frame Transfer (NFT) mode. In this mode of operation, an FT EMCCD acts much like a standard CCD. The capture sequence for this mode is illustrated here:




Figure 23: Capture sequence (NFT mode)

- Step 1** Both **Image** and **Storage** areas of the EMCCD are fully cleared out.
- Step 2** On receipt of a start acquisition command the EMCCD stops the Keep Clean Cycle and an acquisition begins. The image builds up in the **Image** area of the EMCCD until the exposure time has elapsed, at which point the readout process starts.
- Step 3** The first phase of this process is to quickly shift the charge built up in the Image area into the **Storage** area. The time required to move the charge into the Storage area is the same as in Frame Transfer mode
- Step 4** With the image now in the Storage area the captured image is read out. The time taken to read out the image is again the same as in the **Frame Transfer** mode.
- Step 5** On completion of the readout, the EMCCD is completely cleared, ready to acquire the next image. The EMCCD will remain in the Keep Clean Cycle until the end of the accumulation or kinetic cycle time, depending on the acquisition mode, i.e. back to **Step 1**. As at least one Keep Clean Cycle is performed between each exposure, the minimum exposure time is no longer set by the time to read out the image.

---

### 4.4.2.1 Points to note about using an FT EMCCD as a standard EMCCD

- The exposure time, accumulation cycle time and the kinetic cycle time are independent.
- The minimum exposure time is not related to the time taken to read out the image.
- As the captured image is quickly shifted into the Storage area, even in NFT mode, the system can still be used without a mechanical shutter.
- For short exposure times, the image may appear streaked as the time taken to shift the Image area into the Storage area and the exposure time may be of similar magnitude, but much less than a non-frame transfer.
- In conditions where light level falling onto the pixels in the Image area exceeds the pixel well depth of those pixels during the readout of the Storage area image charge blooming can occur vertically along the column contaminating the image being readout.

#### Example:

During a 100  $\mu$ s exposure enough light has fallen on a pixel to register 10000 counts, or 100,000 electrons assuming 10  $e^-$ /count. The image is then shifted into the Storage area. To read out the image, assuming 1,000 x 1,000 pixels, it would take approximately 100 ms at 10 MHz readout rate. This means that during the reading out of the image 10 million counts (10,000 \* 1,000) will have been acquired into the pixel described above. As a pixel saturates at approximately 160,000 electrons this means that the pixel will be over saturated by 60 times. All the excess charge has to go somewhere, and spreads vertically along the EMCCD column. As the clocks in the Image area are not actively shifting the charge, the mobility of the charge will be low, and you may not see any effect. However, when you consider that more than one pixel in any given column could be exposed to 10,000 counts per 100  $\mu$ s, the chance of corrupting data is correspondingly increased. Reducing the amount of light falling on the EMCCD outside of the exposure period and increasing the exposure time accordingly will reduce the possibility of data corruption.

### 4.4.3 Capture Sequence for Fast Kinetics (FK) with an FT EMCCD

Fast Kinetics (FK) is a special readout mode that uses the actual EMCCD as a temporary storage medium and allows an extremely fast sequence of images to be captured. The capture sequence is illustrated here:

- Step 1 Both the Image and Storage areas of the EMCCD are fully cleaned (the Keep Clean Cycle).
- Step 2 The EMCCD stops the Keep Clean Cycle and the acquisition begins. The image builds up on the illuminated sub-area of the EMCCD.
- Step 3 The EMCCD remains in this state until the exposure time has elapsed, at which point the complete EMCCD is clocked vertically by the number of rows for the sub area of the EMCCD.
- Steps 4 & 5 The process is continued until the number of images stored equals the series length set by the user.
- Step 6 At this point the sequence moves into the readout phase by first vertically shifting the first image to the bottom row of the EMCCD. The EMCCD is then read out in the standard method.



---

#### 4.4.3.1 Points to consider when using Fast Kinetics mode

- Light must only be allowed to fall on the specified sub-area. Light falling anywhere else will contaminate the data
- The maximum number of images in the sequence is set by the position of the sub-area, the height of the sub-area and the number of rows in the EMCCD (Image and Storage area)
- There are no Keep Clean Cycles during the acquisition sequence
- Due to the very short exposure times, streaking may be evident

## Section 5: Triggering Information

This section describes the Keep Clean Cycles and Triggering modes for the iXon Ultra and Life 897 models.

### 5.1 Keep Clean Cycles

iXon Ultra and Life cameras have a range of different Keep Clean Cycles that run depending on the actual model and the state the camera is in. The first Keep Clean Cycle runs while the camera is in an idle state, i.e. waiting for the PC to tell it to start an acquisition sequence. The next Keep Clean Cycle runs during an internal trigger kinetics series sequence. The final Keep Clean Cycle runs while the camera is waiting for an external trigger event to occur.

#### 5.1.1 Idle Keep Clean Cycle

When the camera is idle, i.e. not actively acquiring images, it repeatedly runs the **Idle Keep Clean Cycle**. This cycle is composed of a vertical shift, followed by a series of horizontal shifts. The number of horizontal shifts is dependent on the number of elements which make up a row. When the Start command is received from the PC, the camera completes the current "Idle Keep Clean Cycle" and then performs a sufficient number of vertical shifts to ensure both Image and Storage regions are completely charge free, see **Figure 24**. This is enough for any charge in the topmost row to be transferred through the storage area, to the register, and then read out. On completion of this sequence the camera is ready to run the exposure sequence. The exact exposure sequence will depend on several factors including the trigger and the readout modes selected. These will be discussed later in this document.




Figure 24: Idle Keep Clean Cycle - Clocking pattern depends on the acquisition modes and are discussed separately.

## 5.1.2 Internal Keep Clean Cycle

The second type of Keep Clean Cycle is called the **Internal Keep Clean Cycle**. It is performed between individual scans in a kinetic series and is relevant to Non-Frame Transfer Mode combined with either Internal or Software Trigger.

When the user configures a kinetics series acquisition, as well as defining the exposure time and the readout mode, they also define the number of scans to capture and the time between the scans. During the time between individual scans the sensor must be kept free of charge to ensure the data captured is a true reflection of the light that fell on it during the exposure period. The Keep Clean Cycle run during this time is very similar to that described in the Idle Keep Clean Cycle, in that the cycle is one vertical followed by a series of horizontals. In this mode, however, the number of times the cycle is repeated is determined by the cycle time set by the user. The Keep Clean Cycle is completed with a sufficient number of vertical shifts to ensure both the Image and Storage areas are charge free.




Figure 25: Internal Keep Clean Cycle - Once readout of a frame occurs, system repeats a pattern of clocking to keep image and storage regions clear of charge, the number being dependent on the kinetic series cycle time.

### 5.1.3 External Keep Clean Cycle

The third Keep Clean Cycle is the **External Keep Clean Cycle**. This cycle uses a different sequence of horizontal and vertical clocking, as it must be able to respond to external events extremely rapidly, but at the same time keep the image area of the sensor charge free. The External Keep Clean Cycle consists of continuous cycles of one vertical shift, both Image and Storage, followed by reading out one full row. one horizontal shift (see Figure 26). When an external trigger is detected the current cycle completes before the exposure phase starts. It is worth noting that although the External Keep Clean Cycle will complete the current cycle, this will not result in the total loss of signal during this time period, as only one vertical shift will have occurred. For pulsed light of very short time duration, picoseconds (i.e. of the order of one vertical shift), the resultant image may appear to have shifted one row.




Figure 26: The External Trigger Keep Clean Cycle - consists of shifting one row down, and then reading that whole row out. When a trigger event is registered, the current cycle will complete before system goes straight to the exposure phase.

---

## 5.2 Triggering Modes

The iXon Ultra and Life cameras have several different triggering modes. These include **Internal**, **External (and Fast External)**, **External Start**, **External Exposure** and **Software Trigger**. Note also that many of these features require iCam technology within the camera, fuller details of which can be viewed through [andor.oxinst.com](http://andor.oxinst.com)

- In **Internal Trigger** the camera determines the exact time when an exposure happens, based on the acquisition settings entered by the user. This is the most basic trigger mode and requires no external intervention.
- In **External Trigger**, once an acquisition has been started, the camera is placed into "External Keep Clean Cycle", which ensures that charge built up on the CCD is kept to a minimum while waiting for the external trigger event. The External Keep Clean Cycle consists of a continuous sequence of one vertical shift followed by one horizontal shift. Once the External Trigger is received the current Keep Clean Cycle is completed and the exposure phase initiated. The exact nature of the acquisition will depend on the user settings and is explained in more detail in a subsequent section. The external trigger is fed via the Ext Trig input on the camera.
- **Fast External Trigger** is for the most part identical to External Trigger - it differs in only one key aspect. In Fast External Trigger the camera does not wait for a sufficient number of Keep Clean Cycles to have been completed to ensure the image area is completely clean of charge before accepting an external trigger event but, instead, allows a trigger event to immediately start the acquisition process. As a result, Fast External Trigger allows a higher frame rate than standard External Trigger. **NOTE: Once a "sufficient" number of Keep Clean Cycles have been performed, External and Fast External Trigger are identical.**
- **External Start** is a mixture of External and Internal Trigger. In this mode the camera performs a sequence of External Keep Clean Cycles while waiting for one external trigger event to occur before starting the acquisition process. Once this external trigger event has occurred, the camera will switch to internal trigger and the acquisition progresses as if the camera was in Internal Trigger mode.
- **External Exposure Trigger** is a mode of operation where the exposure time is fully controlled by the external trigger input. While the trigger input is high, the CCD accumulates charge in the Image area. When the External Trigger goes low, the accumulated charge is quickly shifted into the Storage area and then read out in the normal manner.
- **Software Trigger** is a mode whereby the camera and software are in a high state of readiness and can react extremely quickly to a trigger event issued via software. This mode is particularly useful when the user needs to control other equipment between each exposure and does not know in advance how long such control will take, or if the time taken changes randomly.

These Triggering modes are explained and illustrated in more detail in the following sections.

## 5.3 Triggering Options in Frame Transfer (FT) Mode

### 5.3.1 Internal Triggering

**Internal Triggering** is the simplest mode of operation, in that the camera determines when the exposure happens. By monitoring the Fire output, the user can determine exactly when the camera is “exposing”.

When the camera is idle, it runs the Idle Keep Clean Cycle (Section 5.1.1). On receipt of the Start command from the PC, the camera completes the current Keep Clean Cycle and then performs sufficient vertical shifts to ensure that the Image and Storage regions are completely free of charge. (i.e. Enough vertical shifts for the top row of the image area, to be brought down through storage area and into the register). The camera then starts its real exposure sequence, for which the timing sequence is illustrated in the Figure 27. At this time, the horizontal clocks are still running, clearing out the register.




Figure 27: Internal Triggering in Frame Transfer Mode - Camera begins in idle keep clean, then once a cycle is interrupted, it finished it before shifting image area into the storage area. At that point the exposure begins. Blue sections: Clocks running continuously. Orange sections: Horz and Vert clocks alternating, running one row down, then reading it out.

The first thing to notice is that the Fire output is high for much of the time. This is because there are no Keep Clean Cycles running between each acquisition, and hence the exposure time starts on completion of the transfer of the Image area into the Storage area. This also has the consequence that the exposure time and the cycle time are closely linked. We have defined the exposure time as the time during which there is no vertical shifting of the image area, which also corresponds to the time during which the Fire output will be high. The other point to note is that the exposure time overlaps the read out of the image.

## 5.3.2 External Triggering

When the camera is idle, it runs the Idle Keep Clean Cycles described in **Section 5.1.1**. On receipt of the Start command from the PC, the camera goes into its External Keep Clean Cycles. The camera will repeat these cycles a minimum of X times, where X is the number of image rows on the sensor, before it will accept any External Trigger events. Once this period is over, the camera continues to run the External Keep Clean Cycles until an External Trigger is received. At that point the current External Keep Clean Cycle is completed, and the camera stops all vertical clocking and waits for the programmed user delay period before starting the read phase. During the readout phase the Image area is transferred rapidly to the Storage area. The Storage area is then read out in the normal way.

Once the readout is complete the camera continues to wait for the next external trigger event. During this period the shift register is continually clocked but the Image and Storage areas are not. On the next trigger the camera again waits for the programmed delay before starting the readout phase. The camera continues in this loop until the number of images requested has been captured. Because the Image area is not cleaned between trigger events, the effective exposure time is the time between events. The User Defined Delay is to allow for the capture of events which occur after the trigger pulse. In the case of the first trigger, the effective exposure time is given by the User Defined Delay since Keep Clean Cycles have been running up until the first trigger. This is in contrast to the subsequent exposure periods which are defined by the time between the external trigger events. Thus, for experimental protocols that involve Continuous Wave (CW) light the first image will be dimmer; some protocols may require that this image is discarded.




Figure 28: External Trigger in Frame transfer Mode

### 5.3.3 iCam Enhanced Trigger Mode

Since all iXon Ultra and Life cameras have iCam technology, the rising edge of the external trigger can occur before the end of the previous read out, provided that the falling edge of the Fire pulse occurs after the readout has completed, i.e. the External Trigger is only accepted up to the 'User Defined Delay Period' before the end of the readout. This enhanced trigger mode will result in a higher frame rate (see Figure 29).




Figure 29: iCam-enhanced' External Trigger in Frame Transfer mode. For illustration only, the external triggers are shown with variable periods. Arm indicates when a trigger could be accepted and goes high based on the user defined delay period. iCam allows a trigger to be accepted while system is still being read out.



## 5.4 Triggering Options in Non-Frame Transfer (NFT) Mode

### 5.4.1 Internal Triggering (NFT)

In **Internal Triggering (NFT)** mode, when the camera is idle, it repeatedly runs the Idle Keep Clean Cycles. When the Start command is received from the PC, the camera completes the current Keep Clean Cycle, and then perform sufficient vertical shifts to ensure the Image and Storage regions are completely free of charge. The camera is then ready to start the real exposure sequence.

The timing sequence is illustrated in **Figure 30**. During exposure, the Fire output is high and there will be no vertical clocking. However, the horizontal register keeps running. On completion of the exposure time the FIRE pulse goes low, and the Image area of the CCD will be shifted into the Storage area. As the acquired signal is now safely placed in the masked off region of the CCD, light still falling on the CCD will not contaminate the acquired image while it is being read out. On completion of the readout the camera will perform the Internal Keep Clean Cycles until the user specified cycle time has elapsed. This process continues until the complete series of acquisitions has taken place.




Figure 30: Internal Trigger in Non-Frame Transfer Mode.

## 5.4.2 External & Fast External (NFT) Triggering

In the **External Triggering** modes, once an acquisition starts, the camera goes into "External Keep Clean" Cycle (see Section 5.1.3). As can be seen from Figure 31, the External Keep Clean Cycle runs continuously until the first external trigger event is detected; at which point the system performs sufficient cycles to ensure the image area is fully cleaned, before it will accept an external trigger, at which point the exposure phase starts. During the exposure there are no vertical clocks running. There will, however, be horizontal clocks to ensure that the shift register continues to be kept clean. Once the exposure time has elapsed the charge built up in the Image area is quickly transferred into the Storage area. From the Storage area the charge is read out as normal. At the completion of the readout the camera restarts the External Keep Clean Cycles.

If the camera is in Fast External Trigger mode it will accept a trigger event immediately and start the next exposure, completing the current keep clean cycle, but not performing the number required to fully clear out the image area. Once this period has passed, Normal and Fast External Triggers operate in the same manner.

Fast External Trigger is useful in those cases where there is very little background light, and the user is looking for the fastest frame rate. With Fast External Trigger, you may see variation in the background contribution to the signal from light that may have been allowed to fall on the sensor during the readout of the previous image. Fast external trigger does not mean that when a trigger is accepted the system will respond quicker than in normal external trigger mode.




Figure 31: External Trigger in Non-Frame Transfer Mode

**NOTE:** It does not matter if the trigger occurs at the early phase of the Keep Clean Cycle, as the light signal will not be lost during the completion of the cycle since only one vertical shift will have occurred. For pulsed light of very short duration (of the order of one vertical shift), the resultant image may appear to have shifted one row.

### 5.4.3 External Exposure (NFT) Triggering

**External Exposure (NFT)** mode is distinct from the triggering modes covered previously, in that the exposure period is fully controlled by the width of the external trigger pulse. The exposure period starts on the positive edge and concludes on the negative edge. The exposure is physically ended by shifting the Image area into the Storage area. The Storage area is then readout in the normal manner.

On completion of the readout, the External Keep Clean Cycle is started again.




Figure 32: External Exposure (NFT) Triggering

## 5.4.4 Software (NFT) Triggering

This mode is particularly useful when the user needs to control other equipment between each exposure and does not know in advance how long such control will take, or if the time taken changes randomly. With Software Trigger, the camera and software are in a high state of readiness and can react extremely quickly to a trigger event issued via software.

In this mode the camera runs the Idle Keep Clean Cycle until the Start command is issued by the PC, which is identical to all the modes previously discussed. On receipt of this command, the camera switches to run the normal Internal Keep Clean Cycle until a Software Trigger command is issued by the PC. This event will start the exposure and readout sequence. On completion of the readout, the camera returns to the Internal Keep Clean Cycle until the next Software Trigger is issued.




Figure 33: Software Triggering in Non-Frame transfer Mode

## 5.5 Triggering Options in Fast Kinetics (FK) Mode

### 5.5.1 Internal (FK)

As **Fast Kinetics** uses both the Image and Storage areas as temporary storage areas, the number of options available is quite limited. The simplest mode is again **Internal Trigger** and, as with the internal trigger modes described previously, the system determines when the acquisition begins and then uses the exposure time defined by the user. On completion of the exposure period the camera performs the number of vertical shifts defined by the user, and then again waits for the exposure period before the next set of vertical shifts.

This process is repeated until the number in the series has been captured, at which point the readout starts. The timing sequence is shown in **Figure 34** and as before the Fire output envelopes the period when no vertical clocking is occurring. You will also see there are no readout cycles or Keep Clean Cycles running during the sequence, hence the very fast kinetic cycle period but limited number of exposures in the series.




Figure 34: Internal Triggering in Fast kinetics Mode.

## 5.5.2 External (FK)

In **External Trigger** mode, a trigger pulse is required to start each scan in the series. The rising edge of the trigger defines the trigger event. The user can delay the start of the vertical shifting relative to the trigger event. After the delay has elapsed, the number of rows (as specified by the user) are vertically shifted. The system then waits for the next trigger to start the next scan. As there is no Keep Clean Cycle running while waiting for the External Trigger, the 'real' exposure time is the time between each trigger. A consequence of this is that, if your experiment has a constant background signal but your trigger period is not fixed, you may see different background levels in your signal. As with internal trigger, the data is only read off the sensor when the capture sequence has completed.




Figure 35: External Triggering in Fast Kinetics Mode

### 5.5.3 External Start (FK)

**External Start** triggering mode is a combination of External and Internal Trigger. At the start of the capture process, the camera runs the External Keep Clean Cycle, waiting for a trigger pulse to be applied to the External Trigger input. On receiving the trigger, the exposure starts. The exposure period is defined by the user. On completion of the exposure period, the camera performs the number of vertical shifts defined by the sub-area height (set by the user) and, then again, waits for the exposure period before the next set of vertical shifts. This process is repeated until the number in the series has been captured at which point the readout starts.




Figure 36: External Start Trigger in Fast Kinetics mode

---

## 5.6 Shuttering (iXon Ultra Only)

The iXon Ultra camera is supplied with a built in bi-stable shutter. This bi-stable shutter requires no power to maintain the open or closed state - so it is well suited to long exposures. The bi-stable shutter also has a longer theoretical life because of the lower energy levels used. The shutter is intended for taking background images and protecting the camera from excessive light and dust. It is not designed to operate at the high frame rates the camera is capable of.

- Under normal operation, the shutter should be set to **Permanently Open** and the shutter open and close times to 0 seconds.
- If you do need the shutter to open and close automatically during your experiment, then set the opening and closing times to 27 mS.
- The maximum continuous operation of the shutter is 2 Hz.

**NOTE:** If the camera is powered down before the software is closed, the shutter could be left in an open position. This is design intent and the shutter will return to the closed position when next commanded to by the software. If the software is exited normally, the shutter will be automatically closed.



## Section 6: Additional Functionality

### 6.1 Count Convert

One of the distinctive features of the iXon Ultra and Life series is the capability to quantitatively capture and present data in units of electrons or photons; the conversion applied either in real time or as a post-conversion step- this is called **Count Convert**. Photons that are incident on pixels of an array detector are captured and converted to electrons. During a given exposure time, the signal in electrons that is collected in each pixel is proportional to the signal intensity. In EMCCDs, the signal in electrons is further multiplied in the EM Gain register. The average multiplication factor is selected in the software from the RealGain™ scale. It can be desirable to directly quantify signal intensity either in terms of electrons per pixel, or in terms of incident photons per pixel. However, during the readout process, array detectors must first convert the signal in electrons (the multiplied signal in the case of EMCCDs) into a voltage, which is then digitized by an Analogue to Digital Converter (ADC). Each Analogue to Digital Unit (ADU) is presented as a 'count' in the signal intensity scale, each count corresponding to an exact number of electrons. Furthermore, the signal value in counts sits on top of an electronic bias offset value. For cameras produced from February 2018 this 'baseline' is clamped at 500 in normal operation modes. Prior to this the value was set at 200 counts and 500 for the EX and EXF sensor models.

Therefore, in order to calculate to the original signal in electrons, the electron to ADU conversion factor must be very accurately stored by the camera (which varies depending on the pre-amplifier gain selection chosen through software). Calculation of the signal as absolute electrons also requires knowledge of the bias offset and the EM Gain. The calculation path is shown in Figure 37 below:




Figure 37: Count Convert calculation path

Furthermore, knowledge of the Quantum Efficiency (QE) at each wavelength, and light throughput properties of the camera window, enables this process to be taken a step further allowing the signal to be estimated in photons incident at each pixel. For this step, the user must input the signal wavelength. In fluorescence microscopy, for example, this would correspond to the central wavelength defined by a narrow band emission filter matched to the fluorophore of interest. If the spectral coverage of the signal on the detector is so broad that the QE curve varies significantly throughout this range, then the accuracy of the incident photon estimation would be compromised.

The Count Convert functionality of the iXon Ultra and Life series provides the flexibility to acquire data in either electrons or incident photons, using both real time and post-process facilities. Since the real time feature is processed in hardware, there is little or no impact on the display rate. With the post-process option, it is possible to record the original data in counts and perform this important conversion to either electrons or photons as a post-conversion step, while retaining the original data.

---

## 6.2 OptAcquire

**OptAcquire** is a unique control interface, whereby a user can conveniently choose from a predetermined list of set-up configurations, each designed to optimize the camera for different experimental acquisition types, thus removing complexity from the extremely adaptable control architecture of the iXon Ultra and Life series. The control architecture is extremely tuneable, meaning the camera can be adapted and optimized for a wide variety of quantitative experimental requirements, ranging from fast single photon counting through to slower scan, 16-bit dynamic range measurements. However, successfully optimizing EMCCD technology is not a trivial exercise, with various set-up parameters directly influencing different camera performance characteristics.

OptAcquire has been designed as a unique interface whereby a user can choose from a predetermined list of eleven camera set-up configurations. A variety of set-up parameters are balanced behind the scenes through the OptAcquire menu. Furthermore, advanced users may wish to create their own additional OptAcquire configuration. Control parameters include:

- **EM Gain** – This parameter has a direct bearing on both sensitivity and dynamic range.
- **Vertical clock speed** – Flexibility in this parameter is critical to optimizing the camera for lowest noise, fastest speed, minimal frame transfer smear or maximum pixel well depth.
- **Vertical Clock Amplitude** – Can be used to compensate charge transfer when the sensor is being 'over-clocked' and also to reduce charge leakage into the image area when there is saturated signal in the frame transfer storage area (e.g. when combining very short exposure with a slow readout speed).
- **Horizontal readout speed** – Ranging between 17 MHz and 0.08 MHz for the iXon Ultra. 17 MHz for faster frame rates, 0.08MHz for best dynamic range. For the iXon Life: 17Mhz and 10Mhz.
- **Pre-amplifier gain** – Trading off reduced digitization noise versus accessing full pixel well depth.
- **EM / Conventional amplifier** – EMCCD operation provides ultimate sensitivity. The iXon Ultra can also run in a traditional high dynamic range CCD mode. This conventional mode is recommended for relatively 'brighter' signals, or when it is possible to apply long exposures to overcome read noise floor.
- **Frame Transfer (overlap)** – Overlapped readout is used to achieve 100% duty cycle, ideal for fastest frame rate measurements without switching exposure time between frames. This mode should be deselected for time-lapse experiments.

## 6.2.1 OptAcquire modes

Pre-defined OptAcquire modes include:

Mode	Description
Sensitivity and Speed (EM Amplifier)	Optimized for capturing weak signal at fast frame rates with single photon sensitivity. Suited to the majority of EMCCD applications.
Dynamic Range and Speed (EM Amplifier)	Configured to deliver optimal dynamic range at moderately fast frame rates. Moderate EM Gain applied.
Fastest Frame Rate (EM amplifier)	For when it's all about speed! Optimized for absolute fastest frame rates of the camera. Especially effective when combined with sub-array/binning selections.
Time Lapse (EM Amplifier)	Configured to capture low-light images with time intervals between exposures. Overlap ('frame transfer') readout is deactivated.
Time Lapse and Short Exposures (EM Amplifier)	Configured to minimize vertical smear when using exposure < 3ms.
EMCCD Highest Dynamic Range (EM amplifier)	Combines EMCCD low-light detection with the absolute highest dynamic range that the camera can deliver. Since this requires slower readout, frame rate is sacrificed.
CCD Lowest Noise / Slow readout (Conventional Amplifier, iXon Ultra models only)	Optimized for slow scan CCD detection with lowest noise floor. Recommended for long exposure applications where slow readout can be tolerated.
CCD Highest Dynamic Range (Conventional Amplifier, iXon Ultra models only)	Optimized for slow scan CCD detection with highest available dynamic range. Recommended for brighter signals OR when it is possible to apply long exposures to overcome noise floor.
CCD noise / readout balance (Conventional Amplifier, iXon Ultra models only)	Optimized for slow scan CCD detection, achieving a balance between noise floor and readout time.
Photon Counting	Configuration recommended for photon counting with individual exposures < 10 sec.
Photon Counting with Long Exposures (> 10 sec)	Configuration recommended for photon counting with individual exposures > 10 sec.

## 6.3 Pushing Frame Rates with Cropped Sensor Modes

The iXon Ultra and Life offer Cropped Sensor Mode, which provides the following advantages:

- Specialized readout mode for achieving very fast frame rates (sub-millisecond exposures) from 'standard' cameras
- Continuous rapid spooling of images/spectra to hard disk
- User selectable cropped sensor size – highly intuitive software definition
- Optically Centred Crop Mode (Live cell super resolution mode) Continuous imaging with fastest possible frame rate from centrally positioned ROIs e.g. for Super Resolution Microscopy
- The iXon Ultra is now available with the complementary OptoMask accessory, which can be used to shield the region of the sensor outside of the cropped area

If an experiment demands fast temporal resolution but cannot be constrained by the maximum storage size of the sensor (as is the case for 'Fast Kinetics Mode' of readout), then it is possible to readout the iXon Ultra in 'Cropped Sensor Mode'. In this mode, the user defines a 'sub-array' size from within the full image sensor area, such that it encompasses the region of the image where change is rapidly occurring (e.g. a 'calcium spark' within a cell). The sensor subsequently "imagines" that it is of this smaller defined array size, achieved through software executing special readout patterns, and reads out at a proportionally faster frame rate. The smaller the defined array size, the faster the frame rate achievable.

In order to use Cropped Sensor modes, one has to ensure that no light is falling on the light sensitive area outside of the defined region. Any light collected outside the cropped area could corrupt the images which were acquired in this mode. For microscopy set-ups this is now aided with an accessory called OptoMask, which is available from Andor.

Cropped Sensor Mode has the end result of achieving a much faster frame rate than that obtainable in a conventional 'sub-array' / ROI readout (during which we would still have to vertically shift the unwanted rows). The frame rate increase is achieved by not reading out (i.e. discarding) the unwanted pixels.




Figure 38: Illustration of Cropped Sensor Mode

The active imaging area of the sensor is defined in such a way that only a small section of the entire chip is used for imaging. The remaining area has to be optically masked to prevent light leakage and charge spill-over that would compromise the signal from the imaging area. By cropping the sensor, one achieves faster frame rates because the temporal resolution will be set by the time it requires to read out the small section of the sensor.

## 6.3.1 Cropped Sensor Mode Frame Rates

In biological imaging, Cropped Sensor Mode can be successfully used to enhance performance, and throughput, in super-resolution 'nanoscopic' applications including STORM and PALMIRA.

Imaging frame rates exceeding 1,000/s can be achieved with a sufficiently small crop area. A series of measurements carried out on the Andor iXon Ultra 897 EMCCD camera, has demonstrated that Cropped Sensor Mode, in conjunction with binning, pushed the speed beyond 4,000 frames per second.

The Imaging frame rate potential of the Andor iXon Ultra 897 EMCCD camera under conditions of Normal and Cropped Sensor Mode readouts are compared in the Tables below:

### Normal Mode

Binning	Array size						
	512 x 512	256 x 256	128 x 128	64 x 64	512 x 100	512 x 32	512 x 1
1 x 1	56	110	212	398	267	708	2,881
2 x 2	109	210	394	699	486	1,141	-
4 x 4	206	385	682	1,109	820	1,615	-

### Cropped Mode (Optically Centred crop mode in brackets)

Binning	Array size						
	256 x 256	128 x 128	64 x 64	32 x 32	512 x 100	512 x 32	512 x 1
1 x 1	111 (174)	595 (569)	1,433 (1,490)	3,533 (3,021)	282	857	11,074
2 x 2	215 (329)	1,094 (1,013)	2,481 (2,325)	5,555 (4,048)	541	1,607	-
4 x 4	405 (593)	1,883 (1,661)	3,906 (3,236)	7,751 (4,878)	1,005	2,865	-

### Note

All measurements are made at 17 MHz pixel readout speed with 0.3  $\mu$ s vertical clock speed. It also assumes internal trigger mode of operation. Frame rates shown are for 'Corner Tethered' ROIs, with 'Optically Centred' ROI frame rates shown within brackets.

EMCCD-based adaptive optics, for which smaller format EMCCD sensors are often used, can benefit from cropped sensor readout. EMCCDs can be flexibly optimized in cropped mode to exceed 2,000 fps. Use of cropped sensor mode opens new possibilities for very fast adaptive optics imaging, enabling the users to reach into several thousands of frames per second.

There is also potential to use cropped EMCCDs for multi-spectral fluorescence confocal scanning, as an alternative to the arrays of PMTs that have traditionally been used in this approach. The > 90% quantum efficiency of the back-illuminated sensor, single photon sensitivity, array architecture and rapid pixel readout speed can be exploited to markedly improve this approach. The laser dwell-time should be set to coincide with the time to expose and readout a short row of approximately 32 pixels - sufficient spectral channels to yield effective un-mixing of several known emitting dyes, resulting in a data cube of 512 x 512 x 32 (spectral), and taking less than 1 second to generate. There is a clear sensitivity advantage of EMCCD pixels over the usually employed PMT-technology, which is ~5-fold in the blue-green and up to tenfold in the red wavelengths.

---

## 6.4 Advanced Photon Counting in EMCCDs (iXon Ultra Models Only)

Photon Counting in EMCCDs is a way to overcome the multiplicative noise associated with the amplification process, thereby increasing the signal to noise ratio by a factor of root 2 (and doubling the effective quantum efficiency of the EMCCD). Only EMCCDs with low noise floor can perform photon counting. The approach can be further enhanced through innovative ways to post process kinetic data. The industry-leading dark current and Clock Induced Charge (CIC) specification of the Andor's back-illuminated iXon Ultra 897 model renders it uniquely suited to imaging by Photon Counting.

Photon Counting can only be successfully carried out with very weak signals because, as the name suggests, it involves counting only single photons per pixel. If more than one photon falls on a pixel during the exposure, an EMCCD (or an ICCD for that matter) cannot distinguish the resulting signal spike from that of a single photon event, and thus the dynamic range of a single frame exposure is restricted to one photon.

**Key Fact – To successfully photon count with EMCCDs, there has to be a significantly higher probability of seeing a 'photon spike' than seeing a dark current / CIC 'noise spike'. The iXon Ultra 897 has the lowest dark current / CIC performance on the market, yielding both lower detection limits and higher contrast images.**

Under such ultra low-light conditions, 'photon counting mode' imaging carries the key benefit that it is a means to circumvent the Multiplicative Noise (also known as 'Noise Factor'). Multiplicative Noise is a by-product of the Electron Multiplication process and affects both EMCCDs and ICCDs. In fact, it has been measured to be significantly higher in ICCDs. The noise factor of EMCCDs is well theorized and measured; to account for it you increase the shot noise of the signal by a factor of square root 2 ( $\sim x1.41$ ). This gives the new 'effective shot noise' that has been corrected for multiplicative noise. The effect of this additional noise source on the overall Signal to Noise ratio can be readily viewed in the S/N plots in the technical note entitled 'EMCCD signal to noise plots'.

Photon Counting Mode does not measure the exact intensity of a single photon spike, it merely registers its presence above a threshold value. It does this for a succession of exposures and combines the individual 'binary' images to create the final image. As such, this mode of operation is not affected by the multiplication noise (which otherwise describes the distribution of multiplication values around the mean multiplication factor chosen). The end result is that low-light images, acquired through this mode of acquisition, are improved by a factor of  $\sim x1.41$  Signal to Noise, compared to a single integrated image with the same overall exposure time.

To successfully photon count with EMCCDs, there has to be a significantly higher probability of seeing a 'photon spike' than of seeing a dark current / CIC 'noise spike'. The lower the contribution of this 'spurious' noise source to a single exposure within the accumulated series, the lower the detection limit of photon counting and the cleaner the overall image will be, as demonstrated in Figure 39 below:






Figure 39: "Photon Counting" vs. Standard "EM-on" imaging for very weak signals

Images A, B and C were recorded under identical illumination conditions, identical exposure times and each with EM Gain set at x1,000. The benefit of photon counting under conditions of low clock induced charge (CIC) is evident.

- Images D and E are derived from a larger number of accumulated images, in order to yield a greater measurable Signal to Noise ratio.
- An identically positioned Region of Interest on each image was used to determine S/N of 3.86 and 6.02 for standard and photon counted images respectively. This factor improvement is in accordance with the theory of Photon Counting circumventing the influence of multiplicative noise (noise factor) in EMCCD signals.

## 6.4.1 Photon Counting by Post-Processing

As a post-processing analysis, the user is able to 'trial and error' photon counting for a pre-recorded kinetic series, trading-off temporal resolution vs SNR by choosing how many images should contribute to each photon counted accumulated image.

For example, a series of 1,000 images could be broken down into groups of 20 photon counted images, yielding 50 time points. If it transpires that better SNR is required, the original dataset could be re-treated using groups of 50 photon counted images, yielding 20 time points.




Figure 40: Schematic illustration of how photon counting can be applied to a kinetic series as a post processing step, affording increased flexibility in 'trial and error' trading temporal resolution vs SNR.




## 6.5 Spurious Noise Filter

In some cases, it may be desirable to filter spurious EM-amplified background events to give as 'black' a background as possible, eradicating any remaining 'salt and pepper' noise. It is important to utilize noise selection and filter algorithms that are intelligent enough to accomplish this task without impacting the integrity of the signal itself.

The Spurious Noise Filter (SNF) function of iXon Ultra offers the user a choice of advanced algorithms to try. SNF can be applied either in real time or as a post-processing step. Like the count convert option, the real time processing of the filters is performed in hardware, thereby providing minimal impact on the display rate. The options available to the user for using Spurious Noise filter are as shown below:

- Median (available in Real Time & Post-Processing)
- Level Above (available in Real Time & Post-Processing)
- Interquartile Range (available in Real Time & Post-Processing)
- Noise Threshold (available in Post-Processing)

These can be selected from the Real Time or Post Process options, e.g.



NOTE: Andor spurious noise filters' options make use of advanced algorithms that offer excellent discrimination of spurious noise events with minimal effect on signal integrity.

---

## Section 7: Troubleshooting

### 7.1 Unit does not switch on

- Check power cord is plugged in and connected correctly to mains supply
- If applicable, replace fuse in the supplied mains cable as detailed in [Section 7.4](#)
- If the unit still does not switch on after the checks above have been carried out, contact Andor Technical Support

### 7.2 Support Device not recognised when plugged into PC

- Choose another USB port
- Check connections

### 7.3 Temperature Trip Alarm sounds (continuous tone)

To protect the camera from overheating, a thermal switch has been attached to the heat sink. If the temperature of the heat sink rises above the predefined limit, the power supply to the cooler will cut off and a buzzer will sound. Should the buzzer sound ensure the following:

#### Air Cooling

- Check that the air vents on the sides of the detector head are not blocked
- There is sufficient clearance (100 mm) around the camera
- The ambient air temperature is not above 30°C
- The fan has not been deactivated (or the speed set too low) in software
- Check that no foreign bodies are obstructing the fan's rotation

#### Water Cooling

- That there is sufficient water flow passing through the camera head

**NOTE: When using water cooling, always use water that is above the dew point of the ambient environment to prevent condensation from occurring.**

The thermal cut-out will not reset until the camera has been powered off and the temperature of the metal-work reaches a predefined limit. Operation of the camera under conditions that cause repeated cut-outs is not recommended, as the thermal switch has a limited number of operations.

---

## 7.4 Camera High FIFO Fill Alarm

It has been observed in some systems that a camera may stop acquiring after approximately 1 - 10 seconds. When this occurs, it has been found to be due to insufficient USB bandwidth. The camera includes a buffer (FIFO) to overcome any short-term bandwidth reductions, however, sustained insufficient bandwidth will always cause the buffer to overflow – regardless of what size of buffer is used. The PC should be able to cope with a sustained USB data transfer of equal to or greater than 60 Megabytes/second. Modern machines should all be able to cope with this.

## 7.5 Other Possible Causes for Stoppages in Acquisition

Two causes of USB 2.0 Cameras stopping acquisition have been found thus far:

1. **Power saving settings in the PC BIOS.** All PC systems are now required to be shipped in power saving mode (under EU legislation). In the Dell T5500 models (and likely all newer systems) there is a setting in the BIOS called "C States Control". Ensure this is disabled, as this saves power by sacrificing USB bandwidth. Other PC manufacturers may have similar settings.
2. **Sharing USB bandwidth with other devices.** An example of this is using USB to RS-232 adapters. These significantly reduce the USB bandwidth available to the point where the camera cannot continuously transmit images. This occurs even if the software that uses these USB to RS-232 devices has not been started.

This problem can only be overcome if the PC being used has dual USB Enhanced Host Controllers. For example, Dell 760 machines have ICH10 family southbridges, which have two "Enhanced Host Controllers" (EHC). Each one of these is able to receive the maximum USB rate from an iXon Ultra (assuming there is sufficient RAM, CPU etc. to process the data).

If the machine has only one EHC it may not be possible to operate a Zyla/Neo/iXon Ultra and any other devices that require significant USB bandwidth.

To ensure the Camera has full USB bandwidth it is best to ensure it does not share its EHC with any other High Bandwidth USB devices i.e. do not connect the USB to RS-232 adapters into the same EHC. These two EHCs can be shared between the front and the back USB ports of the PC and so it may be advantageous to map which physical USB port is associated with which EHC.

# 7.6 Use of Multiple High Speed USB 2.0 I/O on One Camera

On PCs with two or more EHCs it may be beneficial to map each physical USB connector to its EHC. The figure below shows how to do this on a Dell T5500 PC using a (modern) USB memory stick and the application "UVCView.x86.exe" from Microsoft. USB devices can be High Bandwidth (Fastest), Full Bandwidth or Low bandwidth (e.g. mouse, keyboard). A modern USB stick will be High Bandwidth. Run UVCView.exe to monitor the USBs ports usage.

## 7.6.1 Mapping each port

Use a High Bandwidth device such as a modern USB memory stick. High bandwidth devices will only appear in an EHC section, where one can see the USB mass storage device under "Enhanced Host Controller - 3A6A". There is a second EHC called 3A6C. Test each physical port by placing the USB memory stick into it and record which EHC it is connected to. After changing the port, the window should update automatically after several seconds – if it does not, press F5 to refresh it.





Figure 41: Example USB Map Hub and layout



---

## Section 8: Maintenance

### 8.1 Cleaning Camera Exterior

Only use a dry, clean, lint free cloth to clean the painted surfaces. If necessary, use a water diluted mild detergent to lightly dampen the cloth- do not use Isopropyl alcohol, solvents or aerosols.

### 8.2 Regular Checks

The state of the product should be checked regularly, especially the following:

- The integrity of the enclosure
- Any water hoses used
- The AC/DC External Power Supply
- The mains cable



**WARNING: DO NOT USE EQUIPMENT THAT IS DAMAGED.**

#### 8.2.1 Annual Electrical Safety Checks

It is advisable to check the integrity of the insulation and protective earth of the product on an annual basis, e.g. U.K. PAT testing of the PS-90 Power supply.

### 8.3 Replacement Parts

The supplied PS-90 PSU is the only recommended external power supply for use with the iXon cameras. If this unit is faulty or damaged, please contact Andor for a replacement. Depending on the Terms and Conditions of your Warranty, you may be charged for this replacement. There are no user replaceable parts in the camera head- please contact your nearest Andor representative (see [Section 1.1](#)) if required.

### 8.4 Fuse Replacement

The camera itself does not have a fuse. However, if a U.K. (BS 1363) mains lead has been supplied, it contains a fuse, whose characteristics are as follows:

- Rated Current: 5 A
- Rated Voltage: 240 Vac
- Type: BS 1362
- Size: 6.3 × 25.4 mm (¼ × 1 inches) cartridge

For continued protection, always replace with a fuse of the same type and rating.



**WARNING: DO NOT USE EQUIPMENT THAT IS DAMAGED.**

# Appendix A: Technical Specifications

## iXon Ultra and Life 897 System Specifications •<sup>1</sup>

	iXon Life 897		iXon Ultra 897	
Sensor	BV: Back Illuminated, standard AR coated		#BV: Back Illuminated, standard AR coated UVB: Back Illuminated, standard AR with additional lumogen coating EXF: Back illuminated, dual AR coated with fringe suppression NEW #BB: Back-illuminated, blue optimized AR coated	
Active pixels	512 x 512			
Pixel size	16 x 16 µm			
Image area	8.2 x 8.2 mm with 100% fill factor			
Pixel Readout Rate	17 MHz	10 MHz	10 MHz	17 MHz
Minimum temperature, air cooled, ambient 20°C Chiller liquid cooling, coolant @ 10°C, >0.75 l/min	-70°C -80°C	-70°C -80°C	-80°C -100°C	-80°C -100°C
Thermostatic Precision	± 0.01°C			
Triggering	Internal, External, External Start, External Exposure, Software Trigger			
System window type	UV-grade fused silica, Broadband Visible-Near Infrared, 0.5 degree wedge		#BV: UV-grade fused silica, Broadband Visible-Near Infrared, 0.5 degree wedge UVB, EXF: UV-grade fused silica, Broadband Vacuum Ultraviolet-Near Infrared, 0.5 degree wedge #BB: UV-grade fused silica, Broadband Vacuum Ultraviolet-Near Infrared, 0.5 degree wedge	
Blemish specification	Grade 1 sensor from supplier. Camera blemishes as defined by Andor Grade A			
Digitization	16-bit (at all speeds)			
PC Interface	USB 2.0 7			
Lens Mount	C-mount			

## Advanced Performance Specifications •<sup>1</sup>

	Life 897	Ultra 897						
Dark current and background events <sup>•2,3</sup>								
Dark current (e-/pixel/sec) @ -80°C	0.00025	0.00030						
Dark current (e-/pixel/sec) @ max cooling	-	0.00015						
Spurious background (events/pix) @ 1000x gain / -85°C	0.005	0.0018						
Active area pixel well depth	180,000 e <sup>-</sup>							
Gain register pixel well depth <sup>•4,5</sup>	800,000 e <sup>-</sup>							
Pixel readout rates	30, 10 MHz	EM Amplifier: 17, 10, 5 & 1 MHz Conventional Amplifier: 3, 1 & 0.08 MHz						
Read noise (e-) <sup>•6</sup>	<1	EMCCD Amplifier				Conventional Amplifier		
MHz		17	10	5	1	3	1	0.08
Without Electron Multiplication		89	65	37	15	9.6	5.3	2.7
With Electron Multiplication		< 1	< 1	< 1	< 1	-	-	-
Linear absolute Electron Multiplier gain	1 - 1000 times via RealGain™ (calibration stable at all cooling temperatures)							
Linearity <sup>•6</sup>	Better than 99.9%							
Vertical clock speed	0.3 to 3.33 μs (user selectable)							
Timestamp accuracy	10 ns							
NEW SRRF-Stream <sup>+</sup> mode	Optional							

- Figures are typical unless otherwise stated.
- The dark current measurement is averaged over the sensor area excluding any regions of blemishes.
- Using Electron Multiplication the iXon is capable of detecting single photons, therefore the true camera detection limit is set by the number of 'dark' background events. These events consist of both residual thermally generated electrons and Clock Induced Charge (CIC) electrons (also referred to as Spurious Noise), each appearing as random single spikes above the read noise floor. A thresholding scheme is employed to count these single electron events and is quoted as a probability of an event per pixel. Acquisition conditions are full resolution and max frame rate (30 MHz readout; frametransfer mode; 1.1 µs vertical clock speed; x 1000 EM gain; 10 ms exposure; -95°C).
- The EM register on CCD201 sensors has a linear response up to ~400,000 electrons and a full well depth of ~730,000 electrons.
- Readout noise is for the entire system. It is a combination of sensor readout noise and A/D noise. Measurement is for Single Pixel readout with the sensor at a temperature of -75°C and minimum exposure time under dark conditions. Under Electron Multiplying conditions, the effective system readout noise is reduced to sub 1 e<sup>-</sup> levels.
- Linearity is measured from a plot of counts vs. exposure time under constant photon flux up to the saturation point of the system, at 10 MHz readout speed.
- iXon Ultra 888 should work with any modern USB 3.0 enabled PC/laptop, as every USB 3.0 port should have its own host controller. iXon Ultra 888 also ships with a USB 3.0 PCI card as a means to add a USB 3.0 port to an older PC, or as a diagnostic aid to interoperability issues.

# Environmental Specifications

Location to be used	Indoor use only
Altitude limit for air-cooling	Up to 2000 m
Altitude limit for water-cooling	Up to 6000 m
Operating temperature	0°C to +30°C ambient (non-condensing)
Storage temperature	-25°C to +50°C
Operating relative humidity	< 70% non-condensing
Pollution degree	Pollution degree 2. Normally only non-conductive pollution occurs. Occasionally, however, a temporary conductivity caused by condensation must be expected.
Ingress protection rating	IP20
Electromagnetic compatibility	This is a Class A product. In a domestic environment this product may cause electromagnetic interference, in which case the user may be required to take adequate measures to shield adjacent equipment.
Cooling vent clearance	100 mm minimum

# Mechanical Specifications

	iXon 897
Dimensions	See "Appendix B: Mechanical Drawings" on page 82
Weight (camera only*)	3.7 kg [8 lb 3 oz]
Weight (external power supply)	0.65 kg

\*The camera weight is the head only with no cables or pipes attached and without water or coolant..



## Camera Power Specifications

Mains Input for Supplied External Power Supply	100 - 240 V, 43 - 67 Hz
Power Consumption*	Camera Head and External Power Supply (Typ./ Max.): 57.5 W/ 105.8 W Camera Head Only (Typ./ Max.): 46.6 W/ 88.6 W
Voltage Rating	12 V
Current Rating	9 A
Mains Overvoltage category	CAT II. An overvoltage category of CAT II means that the equipment is designed to cope with transient voltages above the rated supply that would be experienced by any product connected to a standard single-phase mains socket in a building.


\*Power consumption based on 888 camera as worst-case

## External Power Supply Requirements

	PS-90
Low Voltage Supply	12 V d.c. $\pm$ 5%
Low Voltage Supply Current	9 A
Low Voltage Supply Cable Plug	Lemo Redel 1P Series 3-pin Female Plug with 170° keying (Part No. PAH.N0.3GL.LC65GZ)
Low Voltage Supply Cable Plug Insertion View	 <p>0 V Return and Protective Earth</p> <p>12 V</p> <p>No connection *</p>
Low Voltage Supply Product Socket	Lemo Redel 1P Series 3-pin Male Socket with 170° keying (Part No. PKH.N0.3GL.AG)
Low Voltage Supply Product Socket Insertion View	 <p>0 V Return and Protective Earth</p> <p>12 V</p> <p>0 V Return and Protective Earth</p>
Ripple	120 mV max.
Safety	Certified to an appropriate IEC standard, e.g. IEC 60950-1, and meet the reinforced insulation from mains requirement of IEC 61010-1
Environmental	Ensure that the EPS meets the environmental specification of the overall product (see above)

\* Ideally also connected to 0 V Return and Protective Earth, but not essential. Not connected on supplied PS-90.

## Appendix B: Mechanical Drawings



Note: iXon Ultra 897 shown. iXon Life model does not have a Camera Link connector.

## Appendix C: Dew Point Graph

The relationship between Relative Humidity and Dew Point at varying Ambient Temperature is shown below. This can be used to calculate the minimum temperature the cooling water should be set to.




Figure 42: Calculation of dew point from ambient temperature and relative humidity

For example, when using an iXon Ultra 897, you will need 10°C cooling water to guarantee performance down to -95°C in the relatively dry atmosphere of an air-conditioned lab, cooling water at 10°C should not present any problems.

However, in humid conditions (such as exist in some parts of the world) condensation may occur, resulting in damage to the head. In such conditions you will have to use warmer water (20°C or even higher if it is very humid). The minimum EMCCD temperature would then be limited to a higher value.

---

## Appendix D: EMCCD Technology

### What is an Electron Multiplying CCD?

Electron Multiplying Charge Coupled Device (EMCCD) is a type of image sensor that is capable of detecting single photon events without an image intensifier (achievable by way of a unique electron multiplying structure built into the chip). Extremely weak signals may be detected above the read noise of the camera at any readout speed. This is important, because the traditional problem of combining sensitivity with speed in standard CCDs is that the two are mutually exclusive, i.e. greater read noise detection limits result from faster pixel readout.

### Does EMCCD Technology Eliminate Read Out Noise?

System noise within modern silicon-based detectors has two primary sources: dark current noise and read noise. The higher the noise floor on a detector the less able it is to read out the extremely weak signals associated with ultra low-light imaging.

With thermoelectric cooling, dark current noise can be reduced to negligible levels. An EMCCD's ability to multiply weak signals above the detector's read noise floor, by applying EM Gain, effectively eliminates read noise at any speed by reducing it to  $< 1 \text{ e-/p/s}$ .

### How Sensitive are EMCCDs?

Two parameters significantly influence detector sensitivity, namely Quantum Efficiency (QE) and system noise. QE is a measure of a camera's ability to capture valuable photons. A high QE results in more photons being converted to photoelectrons within the EMCCD pixels.

Once converted, the photoelectrons in a given pixel must overcome the detection limit or noise floor of the camera, which is set by the system noise. EMCCDs deliver superior sensitivity by maximizing QE and minimizing system noise, through the unique gain control feature. Single photon events are now well within the capabilities of super sensitive EMCCD technology.

### What Applications are EMCCDs suitable for?

EMCCD based detectors have been designed for the most demanding of low-light, dynamic applications. The detection limit is as low as single photons.

These levels of sensitivity are vital for low-light, life science and physical science imaging applications such as single molecule detection, live cell microscopy, weak luminescence detection, or demanding astronomy applications (to name only a few).

### What is Andor Technology's Experience with EMCCDs?

Andor Technology was the first company to introduce an EMCCD based detector in 2000. Since then, the company has led the way in the development of EMCCD detectors, introducing the first back illuminated EMCCD in January 2003.

Andor has continually developed EMCCD detectors and offers the widest range of EMCCD based detectors on the market. Please go to [www.emccd.com](http://www.emccd.com) for further details.

---

## EMCCD Sensor

All EMCCD sensors in the iXon Ultra and Life ranges have a frame transfer architecture. The frame-transfer EMCCD uses a two-part sensor in which one-half of the array is used as a storage region and is protected from light by a light-tight mask. Incoming photons are allowed to fall on the uncovered portion of the array and the accumulated charge is then rapidly shifted into the masked storage region for transfer to the serial output register. While the signal is being integrated on the light-sensitive portion of the sensor, the stored charge is read out. Frame transfer devices have typically faster frame rates than full frame devices, and have the advantage of a high duty cycle i.e. the sensor is always collecting light.

A potential disadvantage of this architecture is the charge smearing during the transfer from the light-sensitive to the masked regions of the EMCCD (although they are significantly better than full frame devices). The smearing is more prevalent when exposure times are closer to the time taken to shift the charge under the mask (in the order of milliseconds). With the iXon EMCCD series, vertical clock speeds can be tuned via the software to deliver the fastest parallel shifts which also results in faster overall frame rates (especially when using sub-array and/or pixel binning readout options).

The EMCCD sensor is capable of detecting single photon events without an image intensifier, achievable by way of a unique electron multiplying structure built into the chip. Traditional CCD cameras offered high sensitivity, with readout noises in single figure  $< 10e^-$  but at the expense of slow readout. Hence, they were often referred to as 'slow scan' cameras. The fundamental constraint came from the CCD charge amplifier. To have high speed operation the bandwidth of the charge amplifier needs to be as wide as possible. However, it is a fundamental principle that the noise scales with the bandwidth of the amplifier, hence higher speed amplifiers have higher noise.

Slow scan CCD's have relatively low bandwidth and hence can only be read out at modest speeds, typically less than 1 MHz. EMCCD cameras avoid this constraint by amplifying the charge signal before the charge amplifier and hence maintain unprecedented sensitivity at high speeds. By amplifying the signal, the readout noise is effectively by-passed and, as such, EMCCD readout noise is no longer a limit on sensitivity (and can often be considered negligible).

Please see **Section 4** for further details on EMCCD technology and sensor architecture.

## Vacuum Housing

Unless protected, cooled CCD sensors will condense moisture, hydrocarbons and other gas contaminants that will attack the CCD surface. If that happens, CCD performance will decline proportionally and will eventually fail. Fortunately, the integrity of the sensor can be preserved by housing it in a protective enclosure. However, it is important to understand that all such environments are not the same and the underlying technology used can seriously impact camera life (and performance).

---

## Outgassing

Outgassing is the release of a gas trapped in a material. It is a problem encountered in high-vacuum applications. Materials not normally considered absorbent can release enough molecules to contaminate the vacuum and cause damage to optical sensors, window coatings, etc.

Even metals and glasses can release gases from cracks or impurities, but sealants, lubricants and adhesives are the most common cause. Left unchecked, cooling performance would steadily degrade and therefore lead to increased dark current. Furthermore, resulting electrochemical reactions would eventually destroy the sensor.

## UltraVac™

A permanent hermetic vacuum head is an essential component of high-end imaging and spectroscopy EMCCD cameras. A permanent vacuum requires not only a hermetic seal, but also low outgassing- which sets the real limit on long-term performance. These criteria are what Andor's UltraVac™ vacuum process uniquely ensures.

Andor's proprietary UltraVac™ process minimizes outgassing, ensuring peak quantum efficiency and cooling will not degrade, even after years of operation. Temperature of the sensor can be reduced significantly (down to -100°C with an enhanced thermoelectric Peltier design) translating into substantially lower darkcurrent and fewer blemishes. This is particularly critical to EMCCD technology, where even a single thermal electron is detected as a spurious noise spike. Elimination of condensation and outgassing means that the system can also use only a single entrance window, with antireflection coating so that QE of the system is maximized. All vacuum processes are carried out in a Class 1,000 clean room. Andor's UltraVac™ is a proven solution with over 10 years of supplying vacuum systems to the field with a negligible failure rate (Mean Time Between Failure (MTBF) of 100 years).

## Thermoelectric Cooler

The iXon Ultra range has a four-stage Peltier cooling assembly, which utilizes the thermoelectric effect to rapidly cool the sensor down to the stable operating temperature. TE coolers have a cold side (in contact with the sensor) and a hot side. Heat must be efficiently dissipated from the TE cooler for effective cooling of the sensor. The iXon Ultra is designed to yield maximum heat dissipation, via either forced air cooling (in-built fan) or water cooling which, in combination with Andor's UltraVac™ vacuum process, results in market-leading cooling performance. A re-circulator or a chiller can be purchased from Andor to provide convenient and effective heat dissipation through water cooling.

The iXon Ultra camera also contains temperature control components, which regulate the cooling of the camera and ensure that a stable temperature is maintained between and throughout measurements.




Figure 43: Ultravac™ metal hermetic vacuum sealing

## Appendix E: Glossary

This glossary provides an overview of the concepts and terminology used in Andor's EMCCD technology.

### Readout Sequence of an EMCCD

In the course of readout, charge is moved vertically into the shift register then horizontally from the shift register into the output node of the amplifier. The simple readout sequence illustrated below (which corresponds to the default setting of the Full Resolution Image binning pattern) allows data to be recorded for each individual element on the EMCCD-chip. Other binning patterns are achieved by summing charge in the shift register and/or the output node prior to readout. For further information on binning, please refer to [Section 4.3.6](#)).




Figure 44: Readout sequence of an EMCCD

1. Exposure to light causes a pattern of charge (an electronic image) to build up on the frame (or Image Area) of the EMCCD-chip.
2. Charge in the frame is shifted vertically by one row, so that the bottom row of charge moves into the shift register.
3. Charge in the shift register is moved horizontally by one pixel, so that charge on the endmost pixel of the shift register is moved into the Gain register.
4. Charge is shifted into the output node of the amplifier.
5. The charge in the output node of the amplifier is passed to the analog-to-digital converter and is read out.
6. Steps 3 and 4 are repeated until the shift register is emptied of charge.
7. The frame is shifted vertically again, so that the next row of charge moves down into the shift register.

The process is repeated from Step 3 until the whole frame is read out.

### Accumulation

**Accumulation** is the process by which data that have been acquired from a number of similar scans are added together in computer memory. This results in improved signal to noise ratio.



---

## Acquisition

An **Acquisition** is taken to be the complete data capture process.

## A/D Conversion

Charge from the EMCCD is initially read as an analogue signal, ranging from zero to the saturation value. **A/D conversion** changes the analogue signal to a binary (digital) number, which can then be manipulated by the computer.

## Background

**Background** is a data acquisition made in darkness. It is made up of fixed pattern noise, and any signal due to dark current.

## Binning

**Binning** is a process that allows charge from two or more pixels to be combined on the EMCCD-chip prior to readout.

## Counts

**Counts** refer to the digitization by the A/D conversion and are the basic unit in which data are displayed and processed. Depending on the particular version of the detection device, one count may, for example, be equated with a charge of 10 photoelectrons on a pixel of the EMCCD.

## Dark Signal

**Dark signal**, a charge usually expressed as a number of electrons, is produced by the flow of dark current during the exposure time. All CCDs produce a dark current, an actual current that is measurable in (typically tenths of) milliamps per pixel. The dark signal adds to your measured signal level and increases the amount of noise in the measured signal. Since the dark signal varies with temperature, it can cause background values to increase over time. It also sets a limit on the useful exposure time. Reducing the temperature of the EMCCD reduces dark signal (typically, for every 7°C that temperature falls, dark signal halves). EMCCD readout noise is low, and in order not to compromise this by shot noise from the dark signal, it is important to reduce the dark signal by cooling the detector. If you are using an exposure time of less than a few seconds, cooling the detector below 0°C will generally remove most of the shot noise caused by dark signal.

## Detection Limit

The **Detection Limit** is a measure of the smallest signal that can be detected in a single readout. The smallest signal is defined as the signal whose level is equal to the noise accompanying that signal, i.e. a Signal to Noise ratio (S/N) of unity. Sources of noise are as follows:

- Shot noise of the signal itself
- Shot noise of any dark signal
- Readout noise

If the signal is small, we can ignore its shot noise. Furthermore, if a suitably low operating temperature and short exposure time can be achieved, the lowest detection limit will equal the readout noise.

---

## Exposure Time

The **Exposure Time** is the period during which the EMCCD collects light prior to readout.

## Frame Transfer

**Frame transfer** is a special acquisition mode that is only available if your system contains a Frame Transfer (FT) CCD or EMCCD. The iXon Ultra and Life models have an FT EMCCD acquisition mode. An FT CCD or EMCCD differs from a standard CCD or EMCCD in 2 ways. Firstly, it contains 2 areas of approximately equal size as shown below:

1. The first area is the **Image Area**, which is located at the top and farthest from the readout register. This is the light sensitive area of the CCD.
2. The second section is the **Storage Area** and is located between the Image Area and the readout register. This section is covered by an opaque mask, usually a metal film, and hence is not sensitive to light.




Figure 45: Frame transfer

The second way in which a FT CCD differs from a standard CCD is that the Image and Storage areas can be shifted independently of each other. These differences allow FT capable CCD or EMCCD devices to be operated in a unique mode where one image can be read out while the next image is being acquired. It also allows them to be used in imaging mode without a shutter. Note: This is only applicable when the camera is running in Accumulate or Kinetic mode.

---

## Noise

### Pixel Noise

The **Pixel Noise** is the variation in the pixel's charge level when exposed to a constant signal flux over a significantly valid period of read levels. The pixel noise is normally expressed as the value of the Root Mean Square (rms) of these variations.

**NOTE: The rms value is approximately x 4 to x 6 smaller than the peak-to-peak variations in the level values read from the pixel.**

Pixel Noise has three main constituents:

- Readout noise
- Shot noise from the dark signal
- Shot noise from the light signal itself

Shot noise cannot be removed due to the laws of Physics. Most simply defined, shot noise is the square root of the signal (or dark signal) measured in electrons.

### Readout Noise

**Readout noise** is due to the amplifier and electronics. It is independent of dark signal and signal levels and is only very slightly dependent on temperature. It is present on every readout, as a result of which it sets a limit on the best achievable noise performance.

### Shot Noise

**Shot Noise** is due to the basic laws of physics and cannot be removed. Any signal, whether it is a dark signal or a light signal, will have shot noise associated with it. Shot noise is a statistical variation in signal level which follows a Poisson distribution. The shot noise relates to the generating signal by the following relationship:-

If the signal or dark signal =  $N$  electrons, then the shot noise is the square root of  $N$ .

You can do nothing about the shot noise of your signal, but by choosing minimum exposures and operating the EMCCD at suitably low temperatures, the dark signal, and consequently the noise from the dark signal, can be reduced.

### Shot Noise from the Signal

**Shot noise from the signal** is caused by dependence on the signal generated by the light falling onto the sensor.

### Shot Noise from the Dark Signal

**Shot noise from the dark signal** is related to the electrons generated within the sensor, Dark Current etc. Therefore, it is dependent on the exposure time and it is very dependent on the temperature.

---

## Calculation of Total Pixel Noise

The **total pixel noise** is not simply the sum of the three main noise components (readout noise, shot noise from the dark signal and shot noise from the signal). Rather, the rms gives a reasonable approximation - thus:

$$\text{total} = \sqrt{\text{readnoise}^2 + \text{darkshot}^2 + \text{sigshot}^2}$$

where:

- total is the pixel noise
- readnoise is the readout noise
- darkshot is the shot noise of the dark signal
- sigshot is the shot noise of the signal

## Fixed Pattern Noise

**Fixed Pattern Noise (FPN)** consists of the differences in count values read out from individual pixels, even if no light is falling on the detector. These differences remain constant from read to read. The differences are due in part to a variation in the dark signal produced by each pixel, and in part to small irregularities that arise during the fabrication of the EMCCD and in part to settling time of the electronics. Since fixed pattern noise is partly due to dark signal, it will change if the temperature changes but, because it is fixed, it can be completely removed from a measurement by background subtraction.

## Quantum Efficiency/Spectral Response

The glossary refers to signals as a number of electrons. Strictly speaking, these are "photoelectrons" created when a photon is absorbed. When a UV or visible photon is absorbed by the detector it can, at best, produce only one photoelectron. Photons of different wavelengths have different probabilities of producing a photoelectron, and this probability is usually expressed as **Quantum Efficiency (QE)** or **Spectral Response**.

QE is a percentage measure of the probability of a single photon producing a photoelectron, while spectral response is the number of electrons that will be produced per unit photon energy. Many factors contribute to the QE of an EMCCD, but the most significant factor is the absorption coefficient of the silicon that serves as the bulk material of the device.

## Readout

**Readout** is the process by which data are taken from the pixels of the EMCCD and stored in computer memory. The pixels, which are arranged in a single row, are read out individually in sequence. Readout involves amplifying the charge on each pixel into a voltage, performing an analogue to digital conversion and then storing the data in computer memory. The time taken to perform this operation is known as the "read time".

## Saturation

**Saturation** is the largest signal the EMCCD can measure. A signal is measured in terms of the amount of charge that has built up in the individual pixels on the EMCCD-chip. A number of factors determine the maximum amount of charge that the EMCCD can handle.

---

## Scans (Keep Clean and Acquired)

The EMCCD is continually being "scanned" to prevent it becoming saturated with dark current (see **Dark Signal**).

- If the scan is being used simply to "clean" the EMCCD (i.e. it is a keep-clean scan), the charge from the EMCCD is discarded
- In an acquired scan, however, the charge undergoes analogue to digital conversion and is acquired into computer memory so that it can be used for subsequent processing and display: it is "read out" (see **Readout**)

Unless the context specifically indicates otherwise, "scan", in this User Guide, generally refers to an acquired scan.

## Shift Register

The **Shift Register** usually consists of a single row of elements (or pixels) running parallel to, and below, the bottom row of light-gathering pixels (the image area) on the EMCCD-chip. The shift register is protected from light by an aluminium mask. The elements in the shift register have a greater capacity to store charge (i.e. a greater "well depth") than the other pixels on the EMCCD-chip.

## Signal To Noise Ratio

The **Signal to Noise Ratio** (commonly abbreviated as **S/N** or **SNR**) is the ratio between a given signal and the noise associated with that signal. Noise has a fixed component and a variable component (shot noise), which is the square root of the signal. Thus, the S/N usually increases (improves) as the signal increases.

The maximum S/N is the ratio between the maximum signal (i.e. the saturation level) and the noise associated with that signal. At near saturation levels the dominant source of noise is the shot noise of the signal.

---

## Appendix F: Other Information

### F.1 Terms and Conditions of Sale and Warranty Information

The terms and conditions of sale, including warranty conditions, will have been made available during the ordering process. The current version for the US is [available here](#), for all other regions (except Japan) please [click here](#).

### F.2 EU/UK REACH Regulation Statement

Andor's EU/UK REACH Regulation statement is available at the [following link](#).

### F.3 Waste Electronic and Electrical Equipment Regulations 2006 (WEEE)

The company's statement on the disposal of WEEE can be found in the Terms and Conditions.



# Appendix G: iXon China RoHS Hazardous Substances Declaration

Name and Content of Hazardous Substances in the Product  
产品中有害物质的名称及含量 产品中有害物质的名称及含量

Hazardous Substance: 有害物质						
Component Name 部件名称	Lead (Pb) 铅	Mercury (Hg) 汞	Cadmium (Cd) 镉	Chromium VI Compounds (Cr <sup>6+</sup> ) 六价铬化合物	Polybrominated Biphenyls (PBB) 多溴化联苯	Diphenyl Ethers (PBDE) 多溴联苯醚
Printed Circuit Board Assemblies (Surface-mount Resistors and Capacitors, and Brass Connectors) 电路板组件 电路板组件 (表面贴装电阻器和电容器 · 以及黄铜连接器)	X	O	O	O	O	O
Hex Stand-offs (see image in table below) 六角隔撑	X	O	O	O	O	O
Screw Locks (see image in table below) 螺丝锁定	X	O	O	O	O	O
All other parts 其余配件	O	O	O	O	O	O

This table was developed according to the provisions of SJ/T 11364  
本表格依据SJ/T 11364 的规定编制  
O - The content of such a hazardous substance in all homogeneous materials of such a component is below the limit required by GB/T 26572  
O - 表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572 规定的限量要求以下  
X - The content of such a hazardous substance in a certain homogeneous material of such a component is above the limit required by GB/T 26572  
X - 表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572 规定的限量要求

This table shows images for parts within the iXon EMCCD Instrument.

