


ANDOR
Andor Marana-X China RoHS Hazardous Substances Declaration

This table is applicable for the following camera models: Marana-4BN1U-SOI Marana-4BN1X-SOI MARANA-4BN1U-SOM MARANA-4BN1X-SOM MARANA-4BN6U-SOI MARANA-4BN6X-SOI MARANA-4BN6U-SOM MARANA-4BN6X-SOM

Component Name 组成部分	Hazardous Substances : 有害物质					
	Lead (Pb) 铅	Mercury (Hg) 汞	Cadmium (Cd) 镉	Chromium VI Compounds (Cr6+) 六价铬化合物	Polybrominated Biphenyls (PBB) 多溴化联苯	Polybrominated Diphenyl Ethers (PBDE): 多溴联苯醚
Sensor 芯片	o	o	o	o	o	o
Printed Circuit Board Assemblies 电路板组件	x	o	o	o	o	o
Fan Assembly 风扇	o	o	o	o	o	o
Cooling Assembly 制冷系统	o	o	o	o	o	o
Heatsink 散热片	o	o	o	o	o	o
Metal Parts 金属部件	o	o	o	o	o	o
Solder 焊接物	o	o	o	o	o	o
Power Supplies 电源	o	o	o	o	o	o
Accessories 配件	x	o	o	o	o	o
Label 标签	o	o	o	o	o	o
<p>o - Indicates that this toxic or hazardous substance contained in all the homogenous materials for this part is below the limit requirement in SJ/T11363-2006. o-表示该部件所有均质材料中所含有害物质均低于SJ/T11363-2006的限量需求。</p> <p>x - Indicates that this toxic or hazardous substance contained in all the homogenous materials for this part is above the limit requirement in SJ/T11363-2006. x-表示该部件所有均质材料中所含有害物质均高于SJ/T11363-2006的限量需求。</p>						
<p>The Environment Protection Use Period is valid only when the product is operated under the conditions defined in the product manual. 环保使用期限仅在产品使用说明书规定的条件下才有效。</p>						

