

TECHNICAL SUPPORT

Operating temperature addendum


OCAM²K & OCAM²S



Table

INTRODUCTION..... 2

OPERATIONAL ENVIRONMENT..... 2

	Name		Date
Prepared by :	S.Lemarchand		10/05/2021
Reviewed by :	P.Feautrier		10/05/2021
Approved by :	JL.Gach		10/05/2021



Introduction

The purpose of this document is to provide the operational environment specification for the OCAM²K and OCAM²S cameras.

Operational environment

Operation conditions (non-condensing condition)	Maximum temperature	50°C
	Minimum temperature	-35°C
Storage conditions	Maximum temperature	80°C
	Minimum temperature	-40°C
Cooling Fluid	Maximum cooling fluid temperature	35°C
	Minimum cooling fluid temperature	Dew point of the room (recommended)*
	Maximum pressure	10 bars



Be careful to the dew point. Even if OCAM² is hermetically sealed and not subject to dew issue, please use cooling fluid at a temperature not below the dew point, otherwise the warranty could be void.



www.first-light-imaging.com

